

Product datasheet for SC211259

TENT5B (NM_052943) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TENT5B (NM_052943) Human 3' UTR Clone
Symbol:	TENT5B
Synonyms:	FAM46B
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_052943
Insert Size:	954 bp
Insert Sequence:	>SC211259 3'UTR clone of NM_052943 The sequence shown below is from the reference sequence of NM_052943. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GCCTATCCACCTGGCTGCCTTGTAACTGACTCAGACCCTGGCCAGAAGGGAAGGGACTGGGCTCACG
GGGTGGGTGGGGCTCCAAGAGTGTGTGGAGGACATGACCAGAGGCGCAGAATGTGCCAGGAGGCCCA
GCACTGCAGGGTTGGGCTTTGATTGAACAGACCAGACTTTCCCGAGCACGAGGCCCTGTGGGCTTAA
TGCCAGCCTGGGCTCTCAGCAGGAGGACCTTTGGTTATTGCATCACACTTTTGGACCAGCATGCCGG
TGGGGTGCAGGCCCATTTGTCTTGGAAGTCAATTTCTTCAGGAGACAAAGCAAGTTTGGGGTGGTAG
CCTTAATGCCCAGCACCTTGATTGATTCCCATGGTCTAAGAGGGTCTTTTGTACTAAGGCCAGCCTGT
GGTGATCCAGCAGCTTTGGAGGGCCTGCAGGCCACACTGACTGATACATGGTATGTTACAGGGACT
TGTCACCTTTGGATCAGGCTGGGCCACTACAGGGCCTGATGGGCCAATTGGGTAATCACTGGGAGGAGG
GAAACTGAGGCCCGAGGCAGATACCTGGGTTTGGGGAGGGCTGGGTTTAGCCCGTTCTCTGAATAACA
TACCCCAAGCTGTGAAGCCTTCTCTTGAGGGTTGAGGTGGGCAGTATACAGGAGCCAAGCAGCCCAAG
GGTGGGTGCAGCCAGGCCCTTCCAGGGTGCCACTCCTTTCCCTTTGTGTTGCAAGGAGCCCTCCT
GCCTGGGAGCTTGTCTTTGGAAGGTGCTGAGCCTGAATTTGCCACAGTCAAGTCAACGTTCCCAACCC
ACTGGCTTGGTTACAGGGCTCTCACATGGGGTAGGGGAGCCATACTCCCAACCCCTTGCGTATTCCTCT
GGGTTTGTGATATTTGCACTCTTAACACCTGCCAATAAAGACGGTCTACACTGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_052943.4</u>
Summary:	Probable nucleotidyltransferase that may act as a non-canonical poly(A) RNA polymerase. [UniProtKB/Swiss-Prot Function]
Locus ID:	115572
MW:	34.3