

Product datasheet for **SC211263**

TMEM203 (NM_053045) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMEM203 (NM_053045) Human 3' UTR Clone
Symbol:	TMEM203
Synonyms:	HBEBP1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_053045
Insert Size:	952 bp
Insert Sequence:	>SC211263 3'UTR clone of NM_053045 The sequence shown below is from the reference sequence of NM_053045. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CTCATGATCCGCGCTGTCGGGTCAACTAGCCTCACCGAGGTGCCGGAGAGGGAGCGCTGGACAACCTAG
AATGTTGACCTCGAGCCGAGGCCCTACTTGCAGCGCACCGGAGGAGAGGCTCTCTAGTCTGAAGGCACC
GCCGGCTTGCGCCGAGCTGAGTGCCGGGTTTCCCTATTCCAATCCTGTTTGAATGGTTTCTTCAGCAG
GGCTTAAAGAGCAGCCTTCATCCTGAAAATGTATTTCTTTTGTAAATGCTTTGAGTAGATAATCCT
GAATTGAGGTCATGAGGAGGCCCCAGGCCAGACAGTCTGAACCCCTCTGACACTTGGAACTGAAT
ATAAGTAAATGTCCAGGTGGAATCTGAGTATTTCTGTGGATCCTGGGAAAGTACTGTTGCACAAAGG
CTGCAAGCTGGACTCAGGAATGTCTCCAACAGCAGCGCTGACCTAAGAGCTCCCTGTGCCGTCTAT
CCAGACCAGACTTCGGTAGATGCCTTTGTAGATCTATCACATGTAAACGAGCTTGTATCTCCTTCCCT
GTGCCACGAGAGAGATTGGCTTTTATTCCAGTCTAGGCAGAGACAGAAGAATGTTGAATAAGAGCAGC
ATTAGAGTCTGTCTGTTATCTGTTGCCAAGAAAAGAACTCTGCTGTCCAGGCACTGCTTGGCTTAC
TATCCAGCAAAGACTGCAGTTTTGTGGACTTTTGACCACCTTGGGCTGGCACTCTTAGCACACCTGAG
ACAGATTTAAGCCTCCCTAAGAGACTGAAGAGAGGAACAGGTGTCAGATACTCATAGGCACTGAGATCT
ACAAATGGGAAGCTTGTGAGTGGCCATCTTTGTTGGCTACGAACCTTTGGTTTGATGCCAGTCAGGTG
CCACATGAGAACCTTTGTGAGATGCAAAATAAGTAAGAGAATGTTTCTGAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_053045.2
Summary:	Involved in the regulation of cellular calcium homeotasis (PubMed:25996873). Required for spermatogenesis (PubMed:25996873).[UniProtKB/Swiss-Prot Function]
Locus ID:	94107
MW:	35.5