

Product datasheet for SC211271

SFRS17A (AKAP17A) (NM_005088) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SFRS17A (AKAP17A) (NM_005088) Human 3' UTR Clone
Symbol:	SFRS17A
Synonyms:	721P; AKAP-17A; CCDC133; CXYorf3; DXYS155E; PRKA17A; SFRS17A; XE7; XE7Y
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_005088
Insert Size:	956 bp
Insert Sequence:	>SC211271 3'UTR clone of NM_005088 The sequence shown below is from the reference sequence of NM_005088. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AGCAGGCACCGCAGTACCTGGAACAGGTAAAGACGGGCACGGCCTCCCCACGGCCTGTCCGGGAAAGAC
CAGGACCTGCTCGAGCCTCCTGGCCGCTCCTTGGCCGCTCTCCGTCCACCCCTGCAAAGCCAAGACCCT
TCTGCAGCCACGAATGTCCACGGAGCCCGCGGAGGAAGGAAGACACCATGCTTTAGAGATCCATCTT
TCTCCACTCACCGCAGCGTACTTGGCACTTCAGTTTCAAACACGTAGTCCTTTAAACTTGATCCGATA
GCTTTAATGCGGCCGGTCTCTCTCAGTCAGGAAAATTGCACAGACCGACAGTCGTGAGGATGGCAGAG
CTGCTGCATTCCCCACACGGGGATTTCTGTGTCTGCTTGGCGACCTCCCTGCGTGACGGCCTAGGAG
GTGCACGGGCCACCATAGTCACACTGGCACTGAAAAGAAAGCGTTGCCCTGGTGATTCTTTCCCCCCCG
TTTGTAATGTAACTGATCAGGAAGTGCAGTTTGGGTGGGATGCCGAATCGTCGTGCTGACATTGAGTC
ACGGATGAGGAAGGTACAAGTCCTTTAAGATCAAACTCAAACGGGCCGTTCTTTCTAAGGTGTCGGTA
TGTGGGGAGTGGTACAAAATGGTCTGATGCTCCTTCAAAACATTCACTTTTTACAACGTCAAGGAATT
AAGCATAAAAAGATTGGTTAAAGCTTTGGTTTCTAGTAAAGGTTAGTGTGTGGTTTTTTTAAAGAA
GCTGTTTTGCTAAATTATTTTACTTGAATGTTTCAAACAGATTTCAAGCTGCAAACTTGTTTTATAA
TCGTTTGTCTTCTCAAGTGAAGCTCAGAAATACCTAAAAATAGCTGTAACGTTTCGCGTTAGGAAAGATG
GTGTTTATTCCAGTTTGATTTTATGGTGAAATAAAATCCTTTTCCAATGAACTAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_005088.3
Summary:	This locus encodes a protein kinase A anchoring protein. The encoded protein is part of the spliceosome complex and is involved in the regulation of alternate splicing in some mRNA precursors. Alternatively spliced transcript variants have been identified for this gene. [provided by RefSeq, Sep 2010]
Locus ID:	8227
MW:	35.7