

Product datasheet for SC211282

TOX2 (NM_001098797) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: TOX2 (NM_001098797) Human 3' UTR Clone
Symbol: TOX2
Synonyms: C20orf100; dj495O3.1; dj1108D11.2; GCX-1; GCX1
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_001098797
Insert Size: 960 bp
Insert Sequence: >SC211282 3'UTR clone of NM_001098797
 The sequence shown below is from the reference sequence of NM_001098797. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CCCAGGGACAAATCGCTCTACCTCACCTTAATCCCCGCCTCCCTACCATCCCTGAGGCTCGCTGGAAGGCA
CTGCTCAGAGCCTGAAGGGCTGACAGCAGAAAAGAGGCCCTGGCCAGAGGCAGGGTGGCCCATCGGAGA
GAGCAGTGACACACCCATTGCCCGGGGGCTGAGTCTCTTCTCAACCTCCACCAGACTCTGCAGAGGC
AGCCCACTGCCACCACCAGCCAAAGAACCTGCAGGAACCTTCCGCCCGCTGACCTGCTTGCTCCAGG
GTAAGTGTGGACCCTGTCTCGCCCTGCGCACGGTACCCTATGTCTGGACACCCGGCCCCAGCTCCAGC
CCCAGCCAGGTGGGCCGCCCTGGCGGGTTCGTTACCAACGGACACCCACCCAGATGCATGGGCCA
GAGGGCCGGCCCCCGGCATAGATGTGCACATCGGTTTTCCAGTGTGAACAAAAGATTACGAAACCTAGA
AACTGTTGGTTCCGTGTAAGTAGTTGACTACGTGTTTTAGAACTGTGCTGAAGACATCTGTAAGACTAT
TTTTATCTTAAGGGAGACGCGCAGAAAAGCGGCTGCCAAACCGTTTCGTATCCTCACAGCAAGGACC
GGACGCTTGCTAGCCACCCGGAGCACTGCTCTCTTTAATCATGTATTATCTATTTTAAATTGCCG
GCGACGACTTTTGTCTATTTATGAAGAAACCTTGAGAACGAAGTTACAGCTTATCCTACCGTGTGTGTG
GTTTTGGGTTTTCGTTTTGGTTTTGGTTTTGAGTCGTTTGCAGCTGTTTCTGGCCCTGGCGAGTGT
CTGTCTTGGTGCCAGTCTTCTCTCAAATCTCTTTATAATAAACTTCTGAAAAGCTGAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001098797.2
Summary:	Putative transcriptional activator involved in the hypothalamo-pituitary-gonadal system. [UniProtKB/Swiss-Prot Function]
Locus ID:	84969
MW:	35.8