

Product datasheet for SC211296

Mitochondrial ribosomal protein L11 (MRPL11) (NM_170738) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Mitochondrial ribosomal protein L11 (MRPL11) (NM_170738) Human 3' UTR Clone
Symbol:	Mitochondrial ribosomal protein L11
Synonyms:	CGI-113; L11MT; MRP-L11
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_170738
Insert Size:	962 bp
Insert Sequence:	>SC211296 3'UTR clone of NM_170738 The sequence shown below is from the reference sequence of NM_170738. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GCTGCCCAAGAAGAAGCTGCCAAGAAGTACCTTGCCCCACCAACTCCCAGATTCAAAGGAGGTAGT
TGCAAAAGCTGTGCCAAGGGGAGGAAGGAGGTACACCAATATGATGATGGTTTTCATGACTTTGAAT
GATATATTTTTGTACATCTAGCTGTATCGAGGCATCAGGCCTGAATAAACATCCTTTCTTACCCACTGT
CATCTTCTATCTTGGGACCTCAGGGAAGGGAACCACTGCTGGTGGTAGCCAATATGGGCAGACATTAGG
AAATTGGGGTCTTCATACTGCTTCTCAGATGTGCAAACTGAGGTTCAAGCAAGCAAAAGTAACTTGCCCA
GGTCCTTTTTTTTTCTTTTTTAAGTTATAATTTGTTTTATTGCCAGGTTCTGAATTGGTGCCAGAGGT
AGGAATAGAACCAGGCCTTCTCATTTGAAATCCTGTGTTCTTCCCTCCCCAGGGTGGGAACTATGGT
TGCTAGGGCTGAGGTGTATTGCTTACTTCTCCATCAAAAATCTCTTTATCTCTGAACAACTTGCCAC
CTCTGCAGAAACAAGTCCTGTGGAGATCTCAGTTCAGCCCCAGGCAGGGTACTACATGGAGTACCCA
GATCCCAGGGGCTACTTCAGGCCCTTCTACTAAGCTCCTGACCTGGCTCTGGGTGAGATGTTGGCTGA
ACATGAACTATTGTAGGGCTGACAGCTACAGACCACAGAAAGAGCAAAAGCTTCAGTAGGTTCTGCCTG
ATGGAGAAGTCTCGGGAAGCCCTTGGCCTTCTATTTCTCATTGTCTTTCAGAGAAGAGAGCAGGAG
TTGGGGGAGAAAGGAGATTGCGATTTACTGACCATCTGTTATGTGCTAGACACTTCATTTTACACTGAG
AAAATATGGAAAGCATAGATAAACACCTAAAAGAAAATAAAAACCTCATGTTATCCCAACACCCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_170738.4</u>
Summary:	This nuclear gene encodes a 39S subunit component of the mitochondrial ribosome. Alternative splicing results in multiple transcript variants. Pseudogenes for this gene are found on chromosomes 5 and 12. [provided by RefSeq, May 2014]
Locus ID:	65003
MW:	36.2