

Product datasheet for **SC211297**

RALYL (NM_001100391) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: RALYL (NM_001100391) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: RALYL
Synonyms: HNRPCL3
ACCN: NM_001100391
Insert Size: 963 bp
Insert Sequence: >SC211297 3'UTR clone of NM_001100391
The sequence shown below is from the reference sequence of NM_001100391. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGTCATGAGCTGTTTCTACAGATAAAGTGA TCTGAAATAACGCATGATGCCACAAGCAGAAAAGAGAA
ACTGTGACAACCCAGAAATGTGAAAGGAGGTTTCTTACTGGACAGCAGCATCTTGGTCAATTTAT
ATAAAAACCCAAAATAAATAAATGGACAGTATTGCTCAGTTTTAGAAAATCCATTTCTTCTATGTTTTA
AGCTGTACAATTGTCAGGTTTTTATGGTTAAATTGTAATGTGTTTTCCCTTTGCTAATTATGTTTT
TTTTTTCAGTCTTAAATGTGAAAGGCATTTATGAATGGTAAGGGAAACACTATATACAAATGTATATT
TGTAAGCTATTTTTATGATTAGCATGTTTCACTGTTGATCATATATAAAGTCAGGTGATATTGCAAT
TCTGTATTTAAAGCTTATTTCCAACAATGTCATGTAAGAAAAGATGCATCTTATGCTAGTTTTTATAAT
TTATTTATAATTTATAGTTTAAAGTACTTCAGATCATAATGATAAAATACTTGAAAAAGTTATTTCT
GCCCTGTATAAGCACCTTTTTATTAATAAAGAATGCAGATATTTAGATGTGATATAATAGTTAAAGA
ACTGTTGGTTTGATCTGTGATTAAGTTGAGCATGCTCCGCTCTACTGAACTAAATGATCCAATTATTAC
TTCAGTCTGGGTATGAGATTCATGGACAAGTAAGGACTAGATTGCCAAGGAAAAGACTGTCTTGCCCT
TGATCCAAAAGTTTAAATTAGTGCATACATCATGTCATTTACCTCCTGTTCTAGGAACTCTCCATT
CCCAAGCATTGCCAGTGTTCAGATAATCTTAGCTGTTGTCTTGTGCTGTGGAATGGAAGAAACCA
TCTTCACAGACCGTAGGAGAATCAACATATAATTTCTTAATAAAATACTGTTTCTTTAAAAACAAA
ACGCGT AAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001100391.3](#)

Locus ID: 138046

MW: 38.4