

Product datasheet for **SC211346**

FGFR1 (NM_023108) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: FGFR1 (NM_023108) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: FGFR1
Synonyms: BFGFR; CD331; CEK; FGFBR; FLG; FLJ99988; FLT2; HBGFR; KAL2; N-SAM; OGD
ACCN: NM_023108
Insert Size: 986 bp
Insert Sequence: >SC211346 3'UTR clone of NM_023108
The sequence shown below is from the reference sequence of NM_023108. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCTGCAGGGATGGGTGGTGGTCTGTATAAGCTTTGAGAGCAACGTCCTACTGGGGCTTTGGGGTCAG
CTACACAAGGAAGGCATTTGGACCCTGCCTTTTCATTGCCGAAACCAGAGCCTTCCACCAAGCGTT
TCCCAGTCTTAGCCCTGTGTCCTGAGTTACGTACGATCTTCTGGCAAATGGGGTGCATGATAAGAGCA
TCTCTTACGAAGAGTTGGAAAAACAAATGCCATATATAAATTCTAAGCCATATGAGGACGAGGAGTAAT
GGCATTTTCTCCTTTTTCTCTCACTCCAGACATTCATTGTCCTGAATGCTCCATTAATCCAGGGA
AGGTAATTGCCTAAATCTCCAGTGGATCTCGCAACAGGAAGGAACCAGAAGCTGGGAAAGTTGTTTACC
TCTTTGTCCCAGAGTTAGACCTCATCCTCCCCTAGCTTAGCTGTCTCAGAGATATACTGGCCCTCCCTT
CTCTTCTTTTGTGCTGGTGGTAAAACGCTCTGTAGGTCATTGGCCACTGTCTCCACTCACACCCC
TGCTCCAGTCTGGAGGGAGTGGGTTAAACACAAATAGAATCCATTTGAAGCAGTGATTCTTTTTT
TTTTTTTTTTTTTTAATCAAATGCTTTGGACTTTTGAAGTCCACTTGTCTGTACTTGTAAAAGGGAA
AGAAGGCCGGGCGCAGTGGCTCACGCCTGTAATCCCAGCACTTTAGATCACTTGAGGTGAGGAGTTGA
GACCAGCCCGGCAACATGGTGAACCCCATCTCTACTAAAAATACAAAATTAGCTGTGCATAGTGGT
TGGCACCTGTAGTCCCAGCTACTCAGGAGGCTGAGGCAAGCTAAGCTTGAACCCAGAAGGCAGAGGT
TGCAGTGAGCTGAGATCACGCCACTGCACTCCAGCCTGGGTGACAGAGTGAGTGTGACTCTCGCTTAAA
AAAAAAAAAAAAAAAAAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_023108.2
Summary:	<p>The protein encoded by this gene is a member of the fibroblast growth factor receptor (FGFR) family, where amino acid sequence is highly conserved between members and throughout evolution. FGFR family members differ from one another in their ligand affinities and tissue distribution. A full-length representative protein consists of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment and a cytoplasmic tyrosine kinase domain. The extracellular portion of the protein interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. This particular family member binds both acidic and basic fibroblast growth factors and is involved in limb induction. Mutations in this gene have been associated with Pfeiffer syndrome, Jackson-Weiss syndrome, Antley-Bixler syndrome, osteoglophonic dysplasia, and autosomal dominant Kallmann syndrome 2. Chromosomal aberrations involving this gene are associated with stem cell myeloproliferative disorder and stem cell leukemia lymphoma syndrome. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. [provided by RefSeq, Jul 2008]</p>
Locus ID:	2260
MW:	36.7