

## Product datasheet for SC211359

### ZNF552 (NM\_024762) Human 3' UTR Clone

#### Product data:

|                           |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | 3' UTR Clones                                                                                                                                                                                                                                      |
| Product Name:             | ZNF552 (NM_024762) Human 3' UTR Clone                                                                                                                                                                                                              |
| Symbol:                   | ZNF552                                                                                                                                                                                                                                             |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                           |
| Vector:                   | pMirTarget (PS100062)                                                                                                                                                                                                                              |
| ACCN:                     | NM_024762                                                                                                                                                                                                                                          |
| Insert Size:              | 988 bp                                                                                                                                                                                                                                             |
| Insert Sequence:          | <p>&gt;SC211359 3'UTR clone of NM_024762</p> <p>The sequence shown below is from the reference sequence of NM_024762. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CAGAGAGTTACAAAAAGAAAGGGCTTATGAGTGCAGTGAGTCCAGTCTCATTAAATATTGGAGAATACA
CACCTGAGTAAGATCCTGTGATTGCAGCAAATGTGGAAATGTTACCCAAAGGTCTGCTCTCCTTGGA
TATTGGAGAGTTCACACTAGAGAAGAGTCCTTACAAGTAAATAAATTTGGCAGTTTTGTAGCCACAC
CTCTGTATTCCTTCAGGATTATAGTTAACTGAATGAAGGCCTTACTAGTGTGACAAATACAGGATAT
TTATCCAGAAGTCTGGCTGCATAACTCACAGGAGAGCTGTCATGTGGGAGGATTGCTACATGGGAGGAT
TTGCTTGAGCCTGGGAGGCAGAGGTTGCAGTGAGCCGAGATTGCAACACTGCAGTCCAGTCTGGACCAC
AGAGTGAGACCCTGTTTCAAAAATAAAATATATAAATGACATCCAATAAAGTATGCATATTTGAGCTGT
AGGATAAGTCTTGAATTTCTTTTTCTTTTTTTTGGAGACAGCCTCGCTCTGTCACCCAGGCTGGAGTAC
AGTGGCATAATCTCGGCTCACTGCAAGCCCTGCCTCCTGGGTTTCATGCCATTCTCCTGCCTCAGCCTCC
CTAATAGCTGGGACTACAGGCAACCACCACCACTCCTACCTTTTTTGTATTTTGTAGAGATGGGGTT
TCACCCTGTTAGCCAGGAAGTTCTCAATCTCCTGACCTCGTGATCTGCCGCCTTGGCCTCTCAAAGTG
TTGGGATTGCAGGCGTGAGCCACTGCTCCTGGCTTGAATTTTCTTATAATCCCATGAAACCATCATCAG
AGTCACAATAATGAATTTATCTAATCTATCACCTTTGTATTACCTGGAGCCTTTATAAAGTTTTTCTC
TATGCAGACAACCATTTACTCTGCTGTATAATACTTTCCATTTTCTAGAATTTTCTGACTGATGT
AATAAAGTGTTTCCTCTTATAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

|                 |                                                                                                                                                                                                                                                                                            |
|-----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components:     | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.                                                                                             |
| RefSeq:         | <u><a href="#">NM_024762.3</a></u>                                                                                                                                                                                                                                                         |
| Summary:        | May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]                                                                                                                                                                                                              |
| Locus ID:       | 79818                                                                                                                                                                                                                                                                                      |
| MW:             | 37.3                                                                                                                                                                                                                                                                                       |