

Product datasheet for **SC211369**

PDHX (NM_001166158) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PDHX (NM_001166158) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PDHX
Synonyms: DLDBP; E3BP; OPDX; PDHXD; PDX1; proX
ACCN: NM_001166158
Insert Size: 986 bp
Insert Sequence: >SC211369 3'UTR clone of NM_001166158
The sequence shown below is from the reference sequence of NM_001166158. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AACCTAGAGAATCCTATCCGACTTGCCTAGTCCTCAAAGATAAGAAGTTGGTGTTCAGCTTAGTTGATT
CAGTAGTTGTTACCAAGAAACATATGTTATAGGAAAACAACCTGGTATTTAAGTATGAAGTGGATGAAA
TGTTTATTTATTTAAGGTGAAAGCATTTGACCCAGGGTGTCTTCATCTCAATTTGGGTTTAAATGTTAT
AGAAATAAATGATGATAAACTCTAACTAAATAAGGAAAGAGAATATTTGGTTACTCAGATCCATTTTAA
ACCTCTGGTGCTGTATAAAGGGAATATTAAGTATGATAAATCAAAGTATATGTTGGCTCATTTGAG
CATTTTGAATATTTGAGAATGTATGATACATGTAATAAATAAATAAATAAATAAATAAATAAATAAATAA
ATGTTGAAGGTAGAAGTATCTCAAAGAGATGGCCATTAACCTAGCAGTGGACCTCACTTTTACAAG
CACTGCTCTAGATATACTTGAAGAATTTAATAGGTACAGAAGTTTATTCTGGATAATAAATAAATAAAGG
ATCACACTGTATTAGGGTTATGGCAACATTATTGAATTTTTTATGTACATAAAGCCATATGTTTAGGG
TGGTTTCTATCTGTCTTGTCTTCACTTATATAAATAAATAAATAAATAAATAAATAAATAAATAAATAA
CATGAATGAAAAGAAATGACATTTCAAAAAAATGGTTCAATGAAAAACTATAGCTAAAATATGTAACCT
TTCTAGGTAAACCGCTTGCCCTTCATCTTGAAGTAAATATATTTAAATAAATTTGTTATCTCTTGCCA
AACATTTTGTAGTGTATTTAAAAACAATAAATAAATAAATAAATAAATAAATAAATAAATAAATAAATAA
TGCAGTAAATGTGAGAAAATATAAACATTTCTATTGATTTTAAATGTTAAAACACAAAAGTTCAATA
AACCTTGAAGTTGATTAGTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001166158.2
Summary:	The pyruvate dehydrogenase (PDH) complex is located in the mitochondrial matrix and catalyzes the conversion of pyruvate to acetyl coenzyme A. The PDH complex thereby links glycolysis to Krebs cycle. The PDH complex contains three catalytic subunits, E1, E2, and E3, two regulatory subunits, E1 kinase and E1 phosphatase, and a non-catalytic subunit, E3 binding protein (E3BP). This gene encodes the E3 binding protein subunit; also known as component X of the pyruvate dehydrogenase complex. This protein tethers E3 dimers to the E2 core of the PDH complex. Defects in this gene are a cause of pyruvate dehydrogenase deficiency which results in neurological dysfunction and lactic acidosis in infancy and early childhood. This protein is also a minor antigen for antimitochondrial antibodies. These autoantibodies are present in nearly 95% of patients with the autoimmune liver disease primary biliary cirrhosis (PBC). In PBC, activated T lymphocytes attack and destroy epithelial cells in the bile duct where this protein is abnormally distributed and overexpressed. PBC eventually leads to cirrhosis and liver failure. Alternative splicing results in multiple transcript variants encoding distinct isoforms.[provided by RefSeq, Oct 2009]
Locus ID:	8050
MW:	38.5