

Product datasheet for **SC211388**

C20orf30 (TMEM230) (NM_014145) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C20orf30 (TMEM230) (NM_014145) Human 3' UTR Clone
Symbol:	C20orf30
Synonyms:	C20orf30; dj1116H23.2.1; HSPC274
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_014145
Insert Size:	971 bp
Insert Sequence:	>SC211388 3'UTR clone of NM_014145 The sequence shown below is from the reference sequence of NM_014145. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TATGATGACATTCAGACTTTGATGACTAGCACCCACCCCATAGCTGAGGAGGAGTCACAGTGGAAGTCTG
TCCAGCTTTAAGATATCTAGCAGAACTATAGCTGAGGACTAAGGAATTCTGCAGCTTGCAGATGTTT
AAGAAAATAATGGCCAGATTTTTGGGTCCTTCCCAAAGATGTTAAGTGAACCTACAGTTAGCTAATTA
GGACAAGCTCTATTTTTCATCCCTGGGCCCTGACAAGTTTTTCCACAGGAATATGTATCATGGAAGAAT
AGAGGTTATTCTGTAATGGAAGTGTTCCTGCCACCACCTCTGTAGAGCTGAGCATTTCTTTTAAA
TAGTCTTCATTGCCAATTTGTTCTGTAGCAAATGGAACAATGTGGTATGGCTAATTTCTTATTATTAA
GTAATTTATTTTAAAAATATCTGAGTATATTATCCTGTACACTTATCCCTACCTCATGTTCCAGTGGA
AGACCTTAGTAAAATCAAAGATCAGTGAGTTCATCTGTAATATTTTTTTTACTTGCTTTCTTACTGACA
GCAACCAGGAATTTTTTATCCTGCAGAGCAAGTTTTCAAATGTAAATACTTCTCTGTTTAACAGTC
CTTGACCATTCGATCCAGTTCACCAAGTAGGTTGGACAGCATATAATTTGCATCATTTTGTCCCTTGT
AAATCAAGATGTTCTGCAGATTATTCCTTAAACGGCCGGACTTTTGGCTGTTTCCTAATGAAACATGTA
GTGGTTATTATTAGAGTTTATAGCCGATTGCTAGCACCTTGTAGTATGTCATCATTTCTGCTCATGAT
TCCAAGGATCAGCCTGGATGCCTAGAGGACTAGATCACCTTAGTTTGATTCTATTTTTAGCTTGCAAA
AAGTGACTTATATTCAAAGAAATTAATGTTGAAATCAAATCCTAGAAATAAATGAGTTAACTTC
AAACA
ACGCGT AAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_014145.5</u>
Summary:	This gene encodes a multi-pass transmembrane protein that belongs to the TMEM134/TMEM230 protein family. The encoded protein localizes to secretory and recycling vesicle in the neuron and may be involved in synaptic vesicles trafficking and recycling. Mutations in this gene may be linked to familial Parkinson's disease. [provided by RefSeq, Mar 2017]
Locus ID:	29058
MW:	37.1