

Product datasheet for SC211398

FBXO3 (NM_012175) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FBXO3 (NM_012175) Human 3' UTR Clone
Symbol:	FBXO3
Synonyms:	FBA; FBX3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_012175
Insert Size:	994 bp
Insert Sequence:	>SC211398 3'UTR clone of NM_012175 The sequence shown below is from the reference sequence of NM_012175. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CGCAGACGCCGTGCTCACGCCTTTTTCAGCAAGCCTTCTGCTGATGGAAGCACTAGGATGATTCTAGG
CTGTAAATAGATTTCTCAATAATGTAAATAACTAAATTGTTCTCTGCATATAGCAGGAAAAGTACAT
GAAATATTGTTTCAGGCCCTGGGTTCTATGTGACACTACATTAGGAATTGGATTGTTGGGTTTGCTTT
GTGTTTTTGAGGTAGAGGAAGAAATGGGAATCTTTTTTTCTCTCCAGGAGTCAGTGGAAGAATAGTT
CTCTAGCTAAGGAATGGACATACCTTTGTTTTAAATATTTTATACTTACAAAAATCTAGAAATGGAGA
GGGAACGTGTTTTGAATAAGGATTTAAATACCTGCACAAGGATAGAGAGAACTATGTGACTCATTCTG
TGAAAAGACTTCTTGCAATTGTGAGTTATTTAGAAATGATCAAAATTTGTAATTAGGCTAATCCATTTA
GTGATTCCTAATATTTGTACTCACAGAGAACTAATTGACTAAACAACCTGAACGCTAGTGGTTTGTC
TTAGACAATCTGTCTTTGAATTTAAAGTCTTTATCGCTAAGACCTTGACTTTAAATTTTTCATCACTAC
AACCTTGAATTTAATTTCAAGTCTTCAACATGATGACCTTGGATTTAATTTAAAGTCTTCAACACTATG
CGCTTTATCATATTATTCACAGATGCATTTTGAATGTAGTATGTAAAAGTATGTAACTGCTGTTTA
TTAACAAAAGATTGTTCAACATCTCATGTAGTTTAAATTTGTAATACTGCTTCTGTTTGTCTCTC
CTTTATACACTTGACTGTCTTTGTGATAAGTGACATGAATTTTATGTTAGGATTAAAGTATGTTTCTG
AAACTTGGATTTTTTTGTAATTATATAATTGAGAGTTAAGAATGAAATCCTTCAAGTGTTAAAACTC
ACATTTTAAAGCAAATTTGGTTCCAA
ACGCGT AAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_012175.4</u>
Summary:	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbxs class. Alternative splicing of this gene generates 2 transcript variants diverging at the 3' end. [provided by RefSeq, Jul 2008]
Locus ID:	26273
MW:	39.6