

Product datasheet for SC211422

HMGXB3 (NM_014983) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	HMGXB3 (NM_014983) Human 3' UTR Clone
Symbol:	HMGXB3
Synonyms:	HMGX3; SMF
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_014983
Insert Size:	966 bp
Insert Sequence:	>SC211422 3'UTR clone of NM_014983 The sequence shown below is from the reference sequence of NM_014983. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GAAGAGGAGGTGGCCGAGTGGCAGAAAGGCCAGGCTGTTGTACAGGGACTACACCATCTCTCAAGCC
ATAGTAAGGCCCTTGCTGAGGCAGAGCTATCCAGGGGACCTGCAGAAGTGGTCTCCTGTGGGGAGGGC
CTCTGACTGCTGGGACTGACCAAGAGCTTCCATTCCCTGAGCATGGTGGGACCCAGGGTCTCAGTTC
TCAACCCTCCAGGGGTGAGGAGTGGTACCAGGAAACCTCTTCTGGCCCCGAGAGAGCACTTGGGGGACA
CGGTATGTTTAATGGAGGGGAGGCTGAGGGAAAGGCTCGTAGCTGGTGGGTTCCTGTGGGCCTGCGGTG
TGGGTGAGGGTGGAGGTCTGGGTGGGCTAAGGCGTGAGCCCCAGCACTAGGTGGGAAGGCTGCTGAGG
TCTCTCCACCCCTGAGGAGCCCTGGTTTCAGCCCCCTCAGTCTGATGAATTGCTTAGCCTGTTGCCTT
TGACTAGGGGCTGGGTGGCCTCATTAACCTAGGGGTCCCTTTGGGCTCTTGATTCTCCCTGAAGGAG
GGATGCATTTCTCTTGTCTTCTGTACCCACATTTGGGGGAAGCTGAGGAGGGAGGAACAGTCAGC
CACAGCTCTCTCCAGCACTGCTCTCCACCCCAAGCTTTGAGGAAGAGCATCCCCTTCTCCTTTCC
CTGGCCACTGCTGCTGCAGCCAATATCCTCTCTGGGCTGGGACCCTCTCCACAGAGGGGATGTGGTCC
CTGGTCATGACATAACCTAGCAGCAGTAGGAAAACTCCCTTCTATGAAGGGGAAGCAGACTGGGCCAT
AAGGAAACAGCAGGACTGGCTCAAGTGCCCAAGGTTTGTAGGCGCTGGGAATTGGCCATGTGTTAAT
TTATTGAGTGGAGTAGGTGGCTTTTTTCCCTCCCTCTCCCCCAACAAGAATAAAGTTTATTAATTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_014983.3
Summary:	This gene is one of the non-canonical high mobility group (HMG) genes. The encoded protein contains an HMG-box domain found in DNA binding proteins such as transcription factors and chromosomal proteins. [provided by RefSeq, Aug 2011]
Locus ID:	22993
MW:	34.1