

Product datasheet for **SC211427**

GPBP1L1 (NM_021639) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GPBP1L1 (NM_021639) Human 3' UTR Clone
Symbol:	GPBP1L1
Synonyms:	SP192
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_021639
Insert Size:	982 bp
Insert Sequence:	>SC211427 3'UTR clone of NM_021639 The sequence shown below is from the reference sequence of NM_021639. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GAAACATCAGATGACGATGCCTGGAAGTAGGCATATAAATGCTCACAGTTAAATCTGACCCAGTAACT
CTGTGTGTTTAGGGAGTATACAAAAGAAATCGTTCTTTTCTTTCTTATGTTGTTGAATACTTCATTC
ACAAGGGAAATAATCATATCCCAAAGAGAGCAATTGGCTTGTGTTTGTGTTTGTATTGTTCTTCCCT
GTTATCTGCTTTATAGAGAGAAGTTTGTGTGGTGGGACAGATTTTTTAAACACACTCACACACACACAC
ACATACACACCCAGTATATATGGGCGATGCACAGGTAGGAGCTGGCAGTGCAGGGAAGAGGAGACT
GGTCTGCAGCAACAGCTTCTACTACCAGCCCTTGGGGCACTCACCCCTGTGATCAAGCAATCATTGTCA
ATGACAAAGTGACTATTGAAGTTATAATTGTATTAATAATGCTAATAATTTGGATATTTATTTTAT
TTTTGGCTGCTCGGGTAACTTTAGCCCTTAACCAAGCATATGTGGGTTTTTTGGTTGTTTTTTTTGT
TTTTTTTTCTTTTCTTTTGGGTACAGCTGTAAATATTTGGATATAGGAAATGTTGTGTTATTCT
TGCAGCCTTGATATTCAGGGTGGATTGTAAATATAAATTTTGTGAGATTTCAAAGATTAAGATTATT
TTGATAACATTATTTACAGATTTAAAAGATGTGGTTATCACAAGTCTCGAGGGGAACTACTGCATAA
AATAACTAACTTGAATAAATATTTTGCATCAGTTTGGATTGTCTGTGTAAGAAGTATATTTCTTTA
AAATTAGATGAGTGGATTCTTTTTTAGGGGAAATATTTTGAAGTTGGGGGAAATGTGAATTCAGAC
ACAAGTGTCTCTATTTTAGCAGATACTTTTTGATGCAGAAAAACCTGTGCAATAAATTACCAGTAA
TAAAGTCAGAAACCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_021639.5</u>
Summary:	Possible transcription factor.[UniProtKB/Swiss-Prot Function]
Locus ID:	60313
MW:	37.6