

Product datasheet for SC211434

ANKRD42 (NM_182603) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ANKRD42 (NM_182603) Human 3' UTR Clone
Symbol:	ANKRD42
Synonyms:	PPP1R79; SARP
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_182603
Insert Size:	1000 bp
Insert Sequence:	>SC211434 3'UTR clone of NM_182603 The sequence shown below is from the reference sequence of NM_182603. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AGTATGAACGACTACGTAGAGAAAAATAGAAATGTCAGCTTGATGAATATCGAGCAGAAGTTGATCAAC
TCAGGGAAACACTGGAAAAAATCAAGTCCAAACTTTGTGGCTATGGTTGGTGTCTTTTAATACAT
TTATTTTCTCAAGAAGTATATACAAGAGTGGCCAAGAGTACAAGCTTTGGGCTGGGCTCTTTGGACT
TGAATCCTGGCTATAGTCTTTTGTATTTGGGCTACTCATTTAACCTTTTGTGCCTCAACTATAAACT
GGAGATGATAACTTATACTTTCTAGAGCTGATGACAGGATTAAAGGAATACACACACACACACACAC
ACACACACACACACACACTCTATATGTAACATAAATTTCTAGATACATACACACATCTGTCTATCT
ATCTTCAAACAGCAGCCTAGTTCATAAGTATTATTGTTAACATTGTACCATAGAGGAGAAAAATGTTTT
CATAAGTAATCAATCAGAATTCTAAGACAGCTAGAGGATATGTATTTTAATATTCTAAATAACCAAA
ATAAAGTGCTTTGAATGGAATCCATATTTCTTTCCATAGGGAAGTTCTTCATCAATCATCATGGATG
AATTAATTCTGTTTGAGGCTTGTCACCTCAAATCTGATCTATTAATATTATAGAATCTTCCTGGCTT
CTTTTATTTTGCACAACTAGTCATTGGTCTCTTTAATAACCACTTTAAAAATATTAATAATAAAAA
CATGTTTGTGTATAGTGTTAAAAACAAGATAGATGTTCTTCTGACACTAAGTCCACTCTCCAGAGG
AAGCTATTGTTAACAATTTAGTACATACCATATTAGACAAAATTCTATTTATCTGCAAAACATGTATGTG
TATTTCTATGCATATGCAAAATAAACCCTTCCTCAGTTTCTCTGGGTTTGGTGTCTCACCAGTGGTT
GATTTTCCAAGAAATATAATTAACATAATAATT
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_182603.4</u>
Locus ID:	338699
MW:	39.5