

Product datasheet for **SC211444**

SLTM (NM_001013843) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLTM (NM_001013843) Human 3' UTR Clone
Symbol:	SLTM
Synonyms:	Met
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001013843
Insert Size:	979 bp
Insert Sequence:	>SC211444 3'UTR clone of NM_001013843 The sequence shown below is from the reference sequence of NM_001013843. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TTTAAGGGTGGACCTCCGCGACGATTCTGAAAAATGAGCTCTCTGCCAAGGTTTTAAGATAATTTATTGA
AATCTCTGTAACTTTACTTGACTACTTATGAAGAGGACCTCTGACTTGCTTGAGAGTTCTGTCAGAC
TTTTCTTTTTAAAAATTAACATGATTGCTTTTCTCAATTTTGGAGAAGATGTTTAAATAGTTCTGTTG
TAACTTTTAATAGTTTTGTGTATCATTCACTTTTTTCTTGCAGCACCGAGGCACATTTGAAAAGATG
GAATTGAAGTCGTTTTGTTAACGCTGTGTGAATATAAAGAGTAGTTTGCAGCTGTGTGGTAGTGGTTT
AATTTGCAGCCTTAGCTCTGTGGTGTCTGGCTCTAGAGTTACTTCTTTTTACCAAGCATTTCAGCCTC
CATTTTGAAGGCTGTCTACACTTAAGAAGTCTTAGCTGTCTAATTTTTAGAGAATAAGATTGTTTCATTG
CATTTCTGAGTATTATGTAACCTATTTTGCAGAAGGTAAGTGTACATTAAAGTGCATCTGTGTATCCTG
GTTTAAAAAATGTAATCTTTTTGAAATAAACCTTCATATTCTGTATAGTTGCTAAAGTGTGAGAAC
CTTTTTAATTGTAATGAGAACCGATTTTCAGTTTAGTGTAGCAGCACACTTGTTCAGGTTTGCATGG
TATGAAACCAAATAGATTCATGAAACCTTGCCATGAGGTTTGTTCACAAGGTTCTTAGACCGAGTTG
TGCAGGTAAGTGCATTTTAGGTAATCTGCACTGTTTGTGATGGATAAATCCATCTCTGGGAATTG
TGTGGGTATTAATGTTTCCATGTTCCCAACTATGTTGAGAAGTGAAAAAACCCAGGTTCTAGATGGG
TGAATCAGTTGGGTTTTGTAAATACTTGTATGTGGGAAGACATTGTTGTCTTTTGTGAAAATAAAAA
TCCACACCTGGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001013843.3</u>
Summary:	When overexpressed, acts as a general inhibitor of transcription that eventually leads to apoptosis.[UniProtKB/Swiss-Prot Function]
Locus ID:	79811
MW:	38.3