

Product datasheet for **SC211447**

L1TD1 (NM_019079) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: L1TD1 (NM_019079) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: L1TD1
Synonyms: ECAT11
ACCN: NM_019079
Insert Size: 986 bp
Insert Sequence: >SC211447 3'UTR clone of NM_019079
The sequence shown below is from the reference sequence of NM_019079. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGAGAATTACTGGGAATAATATACCTTAGCACGCCAGGGTGACTACAAACAATATGCTTTCCTCCCCC
AGCATGCATCCAAAAATCAACAAGTAAACGAAAATACACTTCTACCCAGAAGGATGGACAGCTAATAG
CGTACTTGGGGATGAGGAGCAAGGAATATTACAGATATTACCTAGATGTTAATAAAGGGTATGTTTAAA
AAAAATAGGCTGGTCTCAATGTAGTGAGTTATTTCAAGTTGATCACAGTCAGTTACAGATAGAATTCCTT
GTTTTACTTCCCCCACCACCTCCCTACTGCAGTTGACTAGTCTTTTTAGAAATTTATATTATCTGGCT
GGGCACGGTGGCTCACACCTGTAATCCCAGCACTTTGGGAGGCCGAGGCGGGTGGATCACGAGGTCAGG
AGATCGAGACCATCCTGGCTAACACAGTAAACCCCGTCTACTAAAAATAAAAAAATTAGCCGGGC
ATGGTGGCGGGCCTGTAGTCCCAGCTACTTGGGAGGCTGAGGCAGGAGAATGGCGTGAACCCAGGAG
GCGGAGCTTGAGTGAGCCGAGATCACACCACTGCACTCCAGCCGAGCCACAGAGCGACACTCCGTCT
CAAAAAAAGTATATTGTGCTAGCTTGCTAAGAATAAACTTCTATACTGTTGGGGGAGGGCTGCACCT
GTCAAGATAACCTGTCAATGTAGTAGGAAAACAGGAGGGGACAGTAACAGAAAAGCACGGGAAAAGATG
GCAAGGTTAGTTAAAATAGAAAAGTGCTCAGTTCTCATACCTGTAATCCCAGCAGTTTAGGGGGCCAA
GGCAGGTGGGTCACTTGAGCCAAAGAGTTCAAGGCCATCCTGGGCAATGTGGCGAAAAGTGTCTACAAAA
AAATACAAAAAGAGGAAGAAATGATATTTACAAAGTTGTATCATTTGTGTCATGATGTGACTCGAATATA
TTAAATCTTTTCTTTGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_019079.5</u>
Locus ID:	54596
MW:	36.9