

Product datasheet for SC211449

LRRC17 (NM_005824) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LRRC17 (NM_005824) Human 3' UTR Clone
Symbol:	LRRC17
Synonyms:	P37NB
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_005824
Insert Size:	971 bp
Insert Sequence:	>SC211449 3'UTR clone of NM_005824 The sequence shown below is from the reference sequence of NM_005824. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GAATTCATCGATCCTGGGTCTTTGAGATGAACCCCTGCAAGTAGACTTACGTGAATGATTTTTGCTGTG
CCGCTTTTTTAGGGCTCACACATTTAGAAGAATTAGATTTATCAAACAACAGTCTGCAAACTTTGACT
ATGGCGTATTAGAAGACTTGTATTTTTTGAACCTCTTGCGCTCAGAGATAACCCCTGGAGATGTGACT
ACAACATTCACCTCTACTACTGGTTAAAGCACCCTACAATGTCCATTTAATGGCCTGGAATGCA
AAACGCCTGAAGAATACAAAGGATGGTCTGTGGGAAAATATATTAGAAGTTACTATGAAGAATGCCCA
AAGACAAGTTACCAGCATATCCTGAGTCATTTGACCAAGACACAGAAGATGATGAATGGGAAAAAAC
ATAGAGATCACACCGCAAAGAAGCAAAGCGTAATAATTACTATAGTAGGATAAGGTAGAAATTGTTCTG
ATTGTAATTAGTTTTGTATTTTCTATACTGGTGTAGAAAACATATGTTTACATTTGATTAACCTGTGT
GCCTATTTATGCAGGTAATCCAGCTAAAGGAAGCTTTCTTTAATTATAAGTATTATTGTGACTATTAT
AGTAATCAAGAGAATGCTATCATCCTGCTTGCCTGTCCATTTGTGGAACAGCATCTGGTGATATGCAAT
TCCACACTGGTAACCTGCAGCAGTTGGGTCCTAATGATGGCATTAGACTTTCATAATGTCCTGTATAAA
TGTTTTTACTGCTTTTAGAAAAATAAGAAAAAACTTGGTTCATGTTTACATGCCTTTTCGATAGCTGT
TTGTGCATACTTAAAGATGATCAAAATGATTTTATACAAATGCTGTTATAATAAAATGTCATTCCCTAC
CCCTCTACTTTTTTTCAGTAAGTCATCTTATACATTAATAAATTTCCATTTCTGATTTCTGTGCATTC
CAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_005824.3</u>
Summary:	Involved in bone homeostasis. Acts as a negative regulator of RANKL-induced osteoclast precursor differentiation from bone marrow precursors (By similarity).[UniProtKB/Swiss-Prot Function]
Locus ID:	10234
MW:	37.1