

## Product datasheet for SC211452

### PPIL4 (NM\_139126) Human 3' UTR Clone

#### Product data:

|                           |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | 3' UTR Clones                                                                                                                                                                                                                                      |
| Product Name:             | PPIL4 (NM_139126) Human 3' UTR Clone                                                                                                                                                                                                               |
| Symbol:                   | PPIL4                                                                                                                                                                                                                                              |
| Synonyms:                 | HDCME13P                                                                                                                                                                                                                                           |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                           |
| Vector:                   | pMirTarget (PS100062)                                                                                                                                                                                                                              |
| ACCN:                     | NM_139126                                                                                                                                                                                                                                          |
| Insert Size:              | 988 bp                                                                                                                                                                                                                                             |
| Insert Sequence:          | <p>&gt;SC211452 3'UTR clone of NM_139126</p> <p>The sequence shown below is from the reference sequence of NM_139126. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AAAGATAAAGAAAAATCTAAGTATAGATGAAAGATGAAGAGGCAGAATTGAGAGGCTAACATATTTACT
CTTGTCTAACTTAAGAGTGCCAGGAAAGCAGATGCTTAGATTTTGTGTCAAAGCTTGTTATTTTTTCA
TACTAGGATTATGGTCTTTAGATTAATACTGATTATATAGAGCACGGAAGATAAAGAATTGAACATTT
TCTTTGTACTTTTTTACACTAATTTATTGTTATACATAAATGGTAGTCTTCATTTTTGAAGTCTTA
CATTTTCACTCTTTTTTAAATGAAGTATTCATACTACAAAAATACATAAACGTATATATAAAGGGATA
ATAAATGTAAATATCTGTGTACTCATCAGCCAGCTTAAGATACAGATGTTGTCGACATTTTGAAGTTC
CCTAAGGCCCTCTCCCTCTCAAATAATTATTTGGAATTTTGTGTTTGTCTATTATAGTTTTA
CAACATACGTATGTATCTGTAAGTGAAATGTTAATTTGTATGTTTCTGAATTTTATATAAATGGCAAA
ATGTTTACTTCTGTGACTTTCTTTCATTTTTATTGCTATATAGTATTATATAAATATACTACAACCTAT
TCATTCTTGATGGACATTTGAGTTTTGTTATTTGTTTTGCTGTTGAAACAATGCTGCTGTGAACATT
GTATCTAGCAGTTATATTGCTTGATATATTGATAAAATAATGCCAAATTTGTGTTCTGAAGTGACTGTA
CTAGTTTATACCTCTCTTGCCACTGTGAGTTACCTTTGTTTCATAGCATCTCCAACACTTAAGTTTCAGA
CTTAAAAATTTGCCAACCTAGTGGATGTTAAGTGGTATCTACCTGTGATTTCAGCAAGATTACTGTTA
TCTCTGCTCTTCCATGTGTAAGATGTCATTTTTTCTGGCTGCTTTAAGGTTTTCTGTCACTGGTTG
TAAGAAGTTTGATTATGATGTG
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

|                 |                                                                                                                                                                                                                                                                                              |
|-----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).   |
| Components:     | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.                                                                                               |
| RefSeq:         | <u>NM_139126.4</u>                                                                                                                                                                                                                                                                           |
| Summary:        | This gene is a member of the cyclophilin family of peptidylprolyl isomerases. The cyclophilins are a highly conserved family, members of which play an important role in protein folding, immunosuppression by cyclosporin A, and infection of HIV-1 virions. [provided by RefSeq, Jul 2008] |
| Locus ID:       | 85313                                                                                                                                                                                                                                                                                        |
| MW:             | 38.6                                                                                                                                                                                                                                                                                         |