

Product datasheet for SC211467

DSCC1 (NM_024094) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DSCC1 (NM_024094) Human 3' UTR Clone
Symbol:	DSCC1
Synonyms:	DCC1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024094
Insert Size:	947 bp
Insert Sequence:	>SC211467 3'UTR clone of NM_024094 The sequence shown below is from the reference sequence of NM_024094. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GTTTATAATTCGAGAAGACCCATTTCTTAAAGAACAACGGTCTTTTCTTCAGGACTCAAGTTGCTTTA
TAAAGTTGCTGGATAGAAGAAAAATAACTTGACTTTTATTTCTGACTTTTAGAAACCTAGTCTTCTAA
GGCATTTTTCTCCTCATCTTAAGCATTTTGAACTACTGTTATTTCTTAAATGTTATAGGATGTAGCTT
TAAATTATATTTGACTATGTACACAGTGTAACAAGGACAAAAGAATTCTTCTATATTATGTCTG
AGGTCTCTATGATTTATTTAAAGTTTGGCTTTTATATTTAAATGTCATTTCTTACTTAGGATTTGA
TGTTCTAATCATTAATACACTTAAGTTTCCAGAGAACTAAATAGTAATAGTTTCATAGAAGAGAACTA
TTTATTCACTTAAATATCATTTCTCAAGTGAGTGAGTTCCCTCTACTTTTAGCCTTCCACCCAACTG
GAAGCCTCTAGGTGCTATCAATTATTTATATCCATCGTTTACATCCATGAAATTGGCTGAATAATTACT
CCTCTGCCTGGCGTAGACATGCGCTTTGGGAAAAAACGAGTTTATAATCCTATAATGAAGAATACTGG
CACAGGCAATGCTCACTCGAAAACCTCAAGTAATTTCTAGTTGGTTTTGGAATGCTTGATAAAGTTCTT
TTACAGCTTTATTTCTGATTTGTTTTGGTTTAGATCAAAGTTCAAATTAATTTAACTTAGCTAATG
AACTCATCACCAGGACAGTTGGAGGGGGTAGGCCGAGGTTAAATGGTCCAGCTTTCAAAAATGTTAATG
GCTAATCCATAATTAAAGAAGGTTTAACTGTTACTGAAGTTTACAAGTTTATTGTCATGAACATGAAA
TACAAACACGATGGCTTCGAAATGTCTTTCAATAAATGTTTCTGCATTTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_024094.3
Summary:	CHTF18 (MIM 613201), CHTF8 (MIM 613202), and DSCC1 are components of an alternative replication factor C (RFC) (see MIM 600404) complex that loads PCNA (MIM 176740) onto DNA during S phase of the cell cycle (Merkle et al., 2003 [PubMed 12766176]; Bermudez et al., 2003 [PubMed 12930902]).[supplied by OMIM, Dec 2009]
Locus ID:	79075
MW:	37.2