

## Product datasheet for **SC211496**

### **SYTL2 (NM\_032379) Human 3' UTR Clone**

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** SYTL2 (NM\_032379) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** SYTL2  
**Synonyms:** CHR11SYT; EXO4; SGA72M; SLP2; SLP2A  
**ACCN:** NM\_032379  
**Insert Size:** 1006 bp  
**Insert Sequence:** >SC211496 3'UTR clone of NM\_032379  
The sequence shown below is from the reference sequence of NM\_032379. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
ATGCTTTTGGATTGCCAAGATTTCCAAATGAGCCCAAATCCACTGGCTCCTCCACTGAAAACACTAAA
CCGGTGGAATCTGATCTTGAAAATCTGAGTAGGTGGACAAATCCTCACTTTCTATCTATTGCACCTA
AGGAATACTACACAGCATGTAAAAGTCAATCTGCATGTGCTTCTTTGATTACAAGGCCCAAGGGATTTA
AATATAACAAAATGTGTAATTTGTGACTCTAATATTAATAAGATATTTGAACAAGCTAGGAAAATTGA
ATTTCTGCTGCTGCTTCAAGAAAAAGCTGCCAGAGCATTAAACATGGGGTATTGTTAAGAAGCAAAA
ATGTTCTTGTGGCCATCATGTGTTTCACACCACAATTCTGTGCCACAGTTAAGAGGGTCTGGTACCCT
TGCAGGACCTTTGTAGGTTGTGGAAAAAGTCGCAGAAAGATACTCAAAGTGGAGCAGGGAATGGAGAC
AGACATCAGTGATGATAAAAAAAAAAATGGACCTTAAGAACTATTTACTCTGTAATCTCTAATAAAAT
ATGGAATCCATATTAGGCAATGAGACTGAAACTACTGGTGTCTTTCTGCCTTGAGAAAACAAACAGT
TAAACAAGCCTCAAATGTATTTAGTGCCACCCACTGGCCATAGGTACAATTCAGTTGTTGGCTTGTT
TTGACTTAATTCTAAAATAGGTCTCAAGCTGTATTTTTATGAGTTTATTTTTTAAAACCTGCATAT
ATATGATTGTTTTCTTATAACTTTACTATATGAAAGCAGCATAAGAGTAGTCACAAACATGTTTTGCA
ACAAAGTTTTAATTAGAATGTAAGTTGCTCAGTTATACTGTTCTTCTTATGTATGTAATAATTTTCGTAT
TTTGTA AAAACCTTAGAATAAATTATCATTTGATTTAAATTGTATTAGAAAATTAGCGTGACTTCTCA
TTTTAAATAAAATATTTAGGAATTCTAAACATCTAAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** Sgfl-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_032379.3</a>
<b>Summary:</b>	The protein encoded by this gene is a synaptotagmin-like protein (SLP) that belongs to a C2 domain-containing protein family. The SLP homology domain (SHD) of this protein has been shown to specifically bind the GTP-bound form of Ras-related protein Rab-27A (RAB27A). This protein plays a role in RAB27A-dependent vesicle trafficking and controls melanosome distribution in the cell periphery. Alternative splicing results in multiple transcript variants encoding distinct isoforms.[provided by RefSeq, Jun 2009]
<b>Locus ID:</b>	54843
<b>MW:</b>	39.2