

Product datasheet for SC211520

MBD6 (NM_052897) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MBD6 (NM_052897) Human 3' UTR Clone
Symbol:	MBD6
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_052897
Insert Size:	944 bp
Insert Sequence:	>SC211520 3'UTR clone of NM_052897 The sequence shown below is from the reference sequence of NM_052897. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AAAAACAAGAGGAGGAACTGGCCCCAAGCAGCCATACCTGGAGCTGGATCTGACCCTGATTGGGGAG
AGCTGAGTGCTGAGCCTTGGGAGCCCCGTCAGCCACCTGCACCTGTGGACAGTGGGTGGGGGCACTAC
TCCCCACTCAGAGCACAATGCAACTCCTTCCCCTACAATCCCATCCTGAGCCATTGCAGGGGGCAGGG
AAGTTCACCCCCCCCCACCACCCCCCGCCCCCGAAGCCATGTCAGTAAAAAGCCCTGGGGGGGATG
GTATATGGCCCTTTCCCCACCAGCGCTAAGGGGAACACCCCTTCCCCAGGTCTTTTATTTGTTTAAG
TTATTTTGCACAAATGACTCTTTTATATTTAATTCGATTTTCATTGCCTCCCTTCTTAAAGCCAACAGG
CTCAGTTTACAAACCTGTGAGCTACTGTTGGCTGCTGCCCTCCTTCCCAGTGAAAGGTACAAAGCAATA
AGCATCATGCATCCTCCCTTACCCCTCCAACACCCCTCTGCCTCTGGCTCAGGTTGCTCAAAGCACAG
ATCCTCTCTTACCCCGTCCCCAGGTTTGAACACATAGCCTCATTTCAAGGTGTAGCCAGGTTCCCCCG
ACTTTCCTCTGGGATATAAAAAAGGGGGTAAGGGGGCAAAGAGAGCCCTCTGGGCCTCTCTCCCATAC
ACACTACACTGCCCTTCTCCCCCATCAAAACGCTCAGAGACGTTGTGATGATGCGACTGAGGATTAT
GCAACGTGGTCCAACCGAGCGGCCAGCATGACCAGCTGTCCAGGGGCTGCCTCCTGCCTTTTCTTTTG
TAAAGACAAGACCCTTGGGAGTTTAAATTCTGTTTGTACTTGCCTGTGGGGCCTCCACTGCTTTTCT
ATGGGAGACACTCTTAATTTAACAGATGAGAATATTTTGAACCTCTG
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_052897.4</u>
Summary:	Binds to heterochromatin. Does not interact with either methylated or unmethylated DNA (in vitro).[UniProtKB/Swiss-Prot Function]
Locus ID:	114785
MW:	35.6