

## Product datasheet for **SC211525**

### WDR3 (NM\_006784) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	WDR3 (NM_006784) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	WDR3
Synonyms:	DIP2; UTP12
ACCN:	NM_006784
Insert Size:	2000 bp



[View online »](#)



**Summary:** This gene encodes a nuclear protein containing 10 WD repeats. WD repeats are approximately 30- to 40-amino acid domains containing several conserved residues, which usually include a trp-asp at the C-terminal end. Proteins belonging to the WD repeat family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. [provided by RefSeq, Jul 2008]

**Locus ID:** 10885

**MW:** 76.9