

Product datasheet for **SC211538**

TMEM99 (NM_145274) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMEM99 (NM_145274) Human 3' UTR Clone
Symbol:	TMEM99
Synonyms:	MGC21518
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_145274
Insert Size:	1021 bp
Insert Sequence:	>SC211538 3'UTR clone of NM_145274 The sequence shown below is from the reference sequence of NM_145274. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TGTTTGTTTTCTACGTTTACCTCAACAAGGTACCTTATCATTGTCAGCCTTCATCTCCTGATCCAAA
ATAAAATAAAATGCCACAGGTTACTTGATTTTGCATCTGGAATATGATTATTTAAAGAAATTAATAAT
TAAACTGAAAAATGTAGAGTTTAGAAATCTTTTTAAGGTTTCAAAAAATTTCTAGGAGACAACAAAA
TAGTGTTACTAGATTTAGCAAATAAAATACAGAATGCCAATTAATTTGAATTTCAAATGAACATTGAA
TACTTTTTTTGGTATAACTATATTCATGCAATATTTGGGACATATCTGTTATATTATTCATTGTTTAC
CTGAAATTATAATTTAACTGAGAGTTCTGTATTTTATCTGGCAACCTTATTATGAAATATTTAGGAT
CAATATAAAGGAATAATATTTTTAATTTTTACAACTATTTTAATATAAGTAATTTATTTGATTTTCAT
TGTCATCCTGAGAGTTAAGTTTTGCCTTACAGTCTATGTTACTTTACCTGTGTTTTATGAAAAGACTGT
CTCTTATAATGATAGCCACCACCACCATATCTTGACAAGTATTGACATACTGAACTTTAAATTCATTT
TGTGTAGTTCCCGCCTCTTCATGATTGTTTCTGAGGTAATTCAGCTTCTCCGTCTGGCTTTCCATTTTC
AACTCTAGGTTAGTCTTGTTCAAAGTCAACTCATCTAAGACACTTCGAGGCCCTTGCTATCAGCCTCA
ACACAGTGATAAAGGACAACCTTTACTCTGGCATCTGAACAGAGCCCAATTGGAAGAATCAATACAGCTT
AGATGGAGGTCTCCTCATCCAGTGGTAATAAAACATCCACATAGTGGAGCTTAGAGAGCACAGCACTG
GGATGTGAAGACATGGGTTCAAATTGCACCTCTGCCATGTTGAGCTGTACCTTGGGCAAAATTAATTCT
CTGAACCTTAGTTTTCTCATCTATAAATGGGAATATTGGTTTTTAAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_145274.3</u>
Locus ID:	147184
MW:	40.2