

Product datasheet for **SC211541**

KLC2 (NM_001134775) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: KLC2 (NM_001134775) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: KLC2
ACCN: NM_001134775
Insert Size: 934 bp
Insert Sequence: >SC211541 3'UTR clone of NM_001134775

The sequence shown below is from the reference sequence of NM_001134775. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCTCCCGACGAAGCTCCCTGGTGGGCTAATGTCTGAAGGGGCAGCCAGTCACCAGAGCGCCACCTGGC
ACACCCCTCACCCAGCCCTGCGCATGGGCCTGCTGCTTGCCCGCTGTCTCTCCACAGCCCTG
TCTTTTCTGTTCAATCTCAGGGTAACCTTCTCCCTTGTCATCTCAGCCTGAGCCCTGGAGCTGGGCCT
GCCCACTCCAGCTCCATCCCTTATTTATTCCTTCCAGCAGGGCCCTTCCCTAGGTTCCGGCCAGCAG
GAGGTGCCGGCTGGAGTCTCCACCATAGACTCAGTGGCCTGGCCTCCCCAGACCCAGAGCCAAGAACA
CTAAGCACTCGCCGGCCCTTCGGCACCTCGCCCTCCCTCCCGACTCAACCCGGCCGTTGCTTCTGTAT
ATAGAGAAATAAGTTATTGGCCGCGCCCTCCCTCAGTCCACGGTACTACCCGGCCCTCCCTCGTCC
CTCTTCTAGTGGTACCGCCAGGCCTTAATCACCCCAATCCGTGCGGTGGTATCTCCAGGCTCTACA
TTCTCGGGAGCGCGCCTCCCAAGGGGGTCTGGGACCTTCTCGCGCTCCTCTGGCCTCTGAGGGATG
CGTCTACCCGCGCCATCGCCCGTGGCCAGGACGGGGACCTCCCTTAGTCCGTCTCCACCCGCG
GGCCCTGCCCGCATCCCGGCCTTATGCACTGCCCTCCACCCGCCCCGCCAGGACGCGCCGACCC
CGCCCCGGGCACCGCCACCGAGCCATCTGCCTCGCCTCCCCCACGCCTGCAGTTCTCGCGAGGGG
CGGCGACGGTCCCTGGTGGCAGGAGGGCTCCCTGTTGCGGGTGAGGCGGCTGCTCTCATATTTTC
AGATGTTGCTGTAGAAATAAAGACGGTTTAAATCTGA
ACGCGTAAGCGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001134775.2
Summary:	The protein encoded by this gene is a light chain of kinesin, a molecular motor responsible for moving vesicles and organelles along microtubules. Defects in this gene are a cause of spastic paraplegia, optic atrophy, and neuropathy (SPOAN) syndrome. [provided by RefSeq, Mar 2016]
Locus ID:	64837
MW:	33.1