

Product datasheet for SC211564

NOC3L (NM_022451) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NOC3L (NM_022451) Human 3' UTR Clone
Symbol:	NOC3L
Synonyms:	AD24; C10orf117; FAD24
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_022451
Insert Size:	981 bp
Insert Sequence:	>SC211564 3'UTR clone of NM_022451 The sequence shown below is from the reference sequence of NM_022451. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ACGAAATATTTGAAAACATCACTACACAGTAGAGGAATGAAGTCAGTGGACTTTCTTGATATTTGTG
TGTGCAGATGTACATAAAGATGAGTTGTTAACTTAGGATCTTTCTTTTATACAAGGAAAGCTTCCTA
AGAATGTCTAGGAAGAAGAGGAAGAATGACCTTTGCATGGCACAGGGTCTGCCCTATTCTGAATAT
GTCATTCCATCAAGGAGATCAAAAGCCTTTTTTCTCCCAAGTATTGGAAATTACTTTCTTGATGATG
CTGCCTTTTAAAGCTTACGTACATTATAGTTTTTAAAAAAATCTTTGGACTGGATCTTACTGAAGT
GCAGTTGTCTATATTAATAATAGGCATAGAGCACAGAAAAATCAAGACCATGAGAAGACATTTTACCAT
TTAGCTACTTTTTATACTAAATACTCTTTAAATATTTTATTTCATACTGTGGATGGAATGAGAAG
CATTCTAAATTTGAGTTAATATATTTTATGAAGATATTTGAGAAAAGAAAAAATAGCTTGATTACAG
GTTTCATTGGCTTTTGCTGGATGATCCACCTAAAGAAGTTACCTAATTTGGCCTTTTAAAAAGGTGTTA
GTGTTTATTATAGCTACTTTCAAGGAAAGTTGAATATGATTCTAGTCTCTAAAGTTCTTCACGTTTTT
TGACATTCCCTGGAGGGTGAAGTGGGGAAGAATTGCTCCAGGGTAGAAGAACAGGCCCAAGACTTTACC
ATTCTGATCTAGAGACAAAGGATACTCAATGAGGAGCTTTTTTCCCTCTTGGAACAGGTAATGCTT
TTTCTTATAATAATTATAAAACAGTATTTTATGTAACAGCTATTCCCATATTCTAGGAGTGGCCTA
AGAAATGCGTGTTTCAGTGACTAGATTATAAATATTCTCTATTGTGAATAGTTGAATAAACAGCTGTT
TTTTTCTGCTTCCTA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_022451.11</u>
Summary:	May be required for adipogenesis.[UniProtKB/Swiss-Prot Function]
Locus ID:	64318
MW:	38.5