

Product datasheet for **SC211571**

ST3GAL3 (NM_174967) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ST3GAL3 (NM_174967) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ST3GAL3
Synonyms: DEE15; EIEE15; MRT12; SIAT6; ST3GALII; ST3Gal III; ST3GalIII; ST3N
ACCN: NM_174967
Insert Size: 987 bp
Insert Sequence: >SC211571 3'UTR clone of NM_174967
The sequence shown below is from the reference sequence of NM_174967. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAAGCTGGTAAAAGCTCGCGTCATCATTCTAAGCAGTGGCATCTGAGTGGGCCAGCACATGGCCA
TAGAGGCCAGGCACCACCAGGAGCAGCAGCCAGCACCCTACACAGGAGTCTTCAGACCCAGAGAAG
GACGGTGCCAAGGGCCCCAGGGGAGCAAGGCCTTGGTGGAGCAGCCAGAGCTGTGCCTGCTCAGCAGC
CAGTCTCAGAGACCAGCACTCAGCCTCATTAGCATGGGTCTTGATGCCAGAGGGCCAGCAGGCTCCT
GGCTGTGCCAGCAGGCCAGCATGCAGGTGGTGGGACTGGGAGCAAGGCTGCTGCCGAATCACT
TCTCCAATCAGTGTGGTGTATTATCATTTTGTGAATTTGGGTAGGGGGAGGGTAGGGATAATTTAT
TTTTAAATAAGGTTGGAGATGTCAAGTTGGTTCACTTGCCATGCAGGAAGAGGCCCACTAGAGGGCCC
ATCAGGCAGTGTACCTGTTAGCTCCCTGTGGGGCAGGAGTGCCAGGACCAGCCTGTACCTTGCTGTGG
GGCTACAGGATGGTGGGAGGATCTCAAGCCAGCCCCCTCCAGCTCATGACTGTTTGGCCTTTCTTG
GGGAGAAGGCGGGGTATTCCCCTCACCAGCCCTAGCTGTCCCATGGGGAAACCTGGAGCCATCCCTT
CGGAGCCAACAAGACCGCCCCAGGGCTATAGCAGAAAGAACTTTAAAGCTCAGGAGGGTGACGCCAGC
TCCGCCTGCTGGGAAGAGCTCCCCTCCACAGCTGCAGCTGATCCATAGGACTACCGCAGGCCCGGACTC
ACCAACTTGCCACATGTTCTAGGTTTCAGCAACAAGACTGCCAGGTGGTGGGTTCTGCCTTTAGCCTG
GACCAAAGGGAAGTGAGGCCCAAGGAGCTTACCAAGCTGTGGCAGCCGTCCCAGGCCACCCCATGGA
AGCAATAAAGCTCTTCCCTGA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_174967.4
Summary:	The protein encoded by this gene is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The encoded protein is normally found in the Golgi apparatus but can be proteolytically processed to a soluble form. This protein is a member of glycosyltransferase family 29. Mutations in this gene have been associated with a form of autosomal recessive nonsyndromic cognitive disability as well as infantile epileptic encephalopathy. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2017]
Locus ID:	6487
MW:	35.8