

## Product datasheet for **SC211573**

### ST3GAL3 (NM\_174965) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** ST3GAL3 (NM\_174965) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** ST3GAL3  
**Synonyms:** DEE15; EIEE15; MRT12; SIAT6; ST3GALII; ST3Gal III; ST3GalIII; ST3N  
**ACCN:** NM\_174965  
**Insert Size:** 987 bp  
**Insert Sequence:** >SC211573 3'UTR clone of NM\_174965  
The sequence shown below is from the reference sequence of NM\_174965. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAAGCTGGTAAAAGCTCGCGTCATCACATGATCTAAGCAGTGGCATCTGAGTGGGCCAGCACATGGCCA
TAGAGGCCAGGCACCACCAGGAGCAGCAGCCAGCACCCTACACAGGAGTCTTCAGACCCAGAGAAG
GACGGTGCCAAGGGCCCCAGGGGAGCAAGGCCTTGGTGGAGCAGCCAGAGCTGTGCCTGCTCAGCAGC
CAGTCTCAGAGACCAGCACTCAGCCTCATTGAGCATGGGTCTTGATGCCAGAGGGCCAGCAGGCTCCT
GGCTGTGCCAGCAGGCCAGCATGCAGGTGGTGGGACTGGGAGCAAGGCTGCTGCCGAATCACT
TCTCCAATCAGTGTGGTGTATTATCATTTTGTGAATTTGGGTAGGGGGAGGGTAGGGATAATTTAT
TTTTAAATAAGGTTGGAGATGTCAAGTTGGTTCACTTGCCATGCAGGAAGAGGCCCACTAGAGGGCCC
ATCAGGCAGTGTTACCTGTTAGCTCCCTGTGGGGCAGGAGTGCCAGGACCAGCCTGTACCTTGCTGTGG
GGCTACAGGATGGTGGGAGGATCTCAAGCCAGCCCCCTCCAGCTCATGACTGTTTGGCCTTTCTTG
GGGAGAAGGCGGGGTATTCCCCTCACCAGCCCTAGCTGTCCCATGGGGAAACCTGGAGCCATCCCTT
CGGAGCCAACAAGACCGCCCCAGGGCTATAGCAGAAAGAACTTTAAAGCTCAGGAGGGTGACGCCAGC
TCCGCCTGCTGGGAAGAGCTCCCCCTCCACAGCTGCAGCTGATCCATAGGACTACCGCAGGCCCGGACTC
ACCAACTTGCCACATGTTCTAGGTTTTCAGCAACAAGACTGCCAGGTGGTGGGTTCTGCCTTTAGCCTG
GACCAAGGGAAGTGAGGCCCAAGGAGCTTACCAAGCTGTGGCAGCCGTCCCAGGCCACCCCATGGA
AGCAATAAAGCTCTTCCCTGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_174965.4</a>
<b>Summary:</b>	The protein encoded by this gene is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The encoded protein is normally found in the Golgi apparatus but can be proteolytically processed to a soluble form. This protein is a member of glycosyltransferase family 29. Mutations in this gene have been associated with a form of autosomal recessive nonsyndromic cognitive disability as well as infantile epileptic encephalopathy. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2017]
<b>Locus ID:</b>	6487
<b>MW:</b>	35.8