

Product datasheet for SC211594

NUP160 (NM_015231) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NUP160 (NM_015231) Human 3' UTR Clone
Symbol:	NUP160
Synonyms:	NPHS19
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015231
Insert Size:	1017 bp
Insert Sequence:	>SC211594 3'UTR clone of NM_015231 The sequence shown below is from the reference sequence of NM_015231. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CGGGATTATTATATCGTCGGACCTTGTATTTGGATTGTCACCTAGCCTTTGTAACCGCTTGGTGCCT
CTTAGGACTTAAGACTACCCTACAGGAACCCTGTACTCAAGGCCGATTTTGTAACTGTAAATGATGTG
TACAACATTCAAGTCTGCATTCTGCACAAGATAGGAGGGCGGAAGAGTCAGAGGACCCTGTGCTTGCTG
GTGGTGCTAACACAATTTCTGGTGTTCAACCTTGGTCTCAAATAGCTGCTTTTGTATATGATTACAGAG
CTTTTTTAGAGTTTATATTTTTTAAACTACCGAAGACATTCATTATCTGCAAATTAAGACTCACCTTC
ACTTTCCAAAATAGCTGAGGGTTGTTGGCTTGTGTAGCTGACCACCAAAAGCAGTCACTGCAAATCTT
TTAATTCCTCCCTATCACCTTTTGTATTTTAATGCAATTATTTTGGTCCAGAAGTACCTGTATTTCT
GTATTGTACACAAAAGCTAATAATTTTGTGACTTTTTATTTATTTTGGAGGTTTTATATGATCTTCAA
TTGAGTATTAAATAATTTGCCTAGATTAAAGCTAAATGATGACCAGCTAATTAAGAAGATATTTTGA
ATCTGGTTCTGAGCTAAAGTTGAGTAAATTCCTAGCTAAGAAAAAATTGGAAATCCATCATCTATATTA
GCAACAGATTCTCAGAGTAAATTTGTTAACTTCTATGATTTATGATAATCAAGCTGGACTTGATCATACA
AGTTAGTCTCATAATGTATTGGACCAAAATGTTAACTTCATTGGTCAGATTTAGAAGCATTATGCTCA
CAAGTTTTGGGAAAGTGAAAAATAATAAATCATCTTGGATTTTATTCTGTATATTAATAATTTATCTTT
TAAGGAAACAATCTGTATACTACTTGCTGTATAGCCTTTTGACCCTTCTTGAGTTTTTCAGAAGCCTT
TAATTTTATACTTTCAATACCATATTTACATTATATACTTTAATTAACAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_015231.3</u>
Summary:	NUP160 is 1 of up to 60 proteins that make up the 120-MD nuclear pore complex, which mediates nucleoplasmic transport.[supplied by OMIM, Apr 2004]
Locus ID:	23279
MW:	39.4