

Product datasheet for **SC211613**

GCIP interacting protein p29 (SYF2) (NM_015484) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GCIP interacting protein p29 (SYF2) (NM_015484) Human 3' UTR Clone
Symbol:	GCIP interacting protein p29
Synonyms:	CBPIN; fSAP29; NTC31; P29
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015484
Insert Size:	1020 bp
Insert Sequence:	>SC211613 3'UTR clone of NM_015484 The sequence shown below is from the reference sequence of NM_015484. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGAATTTGAAAGAGGAACAGCTGTCTAATCCCTTCAAGAACTGTTTATAGAAGCTTGAGAATGGGGT
AAAAATTTCTGCTAGCAAAATCAAGTTCTTTTGGAAATTTATCAGTAATCCAGAATTTAGTAGCCAT
GCCTTCTCACTCAGCATTAGAAATAAAAATGTGGTTTCTTAAACGTATATCCTTTCATGTATATTTCC
ACATTTTGTGCTTGGATATAAGATGTATTTCTGTAGTGAAGTTGTTTTGTAATCTACTTTGTATACA
TTCTAATTATATTATTTTCTATGTATTTAAATGTATATGGCTGTTAACTTTGAAGCATTTTGGGC
TTAAGATTGCCAGCAGCACATCAGATGCAGTCATTGTTGCTATCAGTGTGGAATTTGATAGAGTCTA
GACTCGGGCCACTTGGAGTTGTGTACTCCAAGCTAAGGACAGTGATGAGGAAGATGGCAGTGGCCACC
GGAGGACTGGAGCAGTCCCTCCTCATGGCGGCCTGTGACCAAGGTCGGGGAGGAGTGGAGCTATCCTTC
CATGATCTGATCATGTACAGTCCCTTTTTAAAAAGCAATAAATGCTTGGGATTAGAATTTCTAACATT
TTTAGTCTATCATTTAAAACAAACATTTGTATATACATTTGTTTATATACATACATGCTTCCATTCT
ACAAACACTGAAAGAATATTTTATTATGGGAAATTCATACACCCAATACCCCTGGTGAGTCCCATATA
CCCAGTACCCTGTTCAACAATCATTAACACTCATGGCCAGTCTGCTTCATCTTTATTCCATTTACTTTC
TTCATGACTCTTATTTTGGACGCATATACTAGATATATAATTTTCATCTATAATTATTTTAGTATGTATCT
CTGAAAGATAAAGACTTAAATGACCACAATACTATTATCACTGCAATAGATTTATTTTTATAAAACTCC
AGAGGGAAAGGTTGTTATTTTTCTCTATGTTGAATTAATGCTAAATATCTAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG

```

Restriction Sites: Sgfl-Mlul



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_015484.5
Summary:	This gene encodes a nuclear protein that interacts with cyclin D-type binding-protein 1, which is thought to be a cell cycle regulator at the G1/S transition. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]
Locus ID:	25949
MW:	38.9