

Product datasheet for **SC211628**

ZNF596 (NM_001042415) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ZNF596 (NM_001042415) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ZNF596
ACCN: NM_001042415
Insert Size: 1007 bp
Insert Sequence: >SC211628 3'UTR clone of NM_001042415

The sequence shown below is from the reference sequence of NM_001042415. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ACTCACACTAAAAAGCAATGAATATGTAGAATCATCAGCTGTAGCGTTAACTAAATACACCAAGG
ACAAACATACTACAGGAATATTATGTCTGTAATCAGTGTGAAAAGCCTTTATTTATATTTACCCTTT
GCTCAACCTAAATGAATCAAGGTAGAGAGAATCCAGATGTATTTAATGTTTATGGCACAACCTTCAAG
CTCTAGGCTGACCATATACAACGTGAGAGAATGAAACTATAGATCAAAGGAATGTGGAGGAGTCTTCAT
CCACAGCTCTGTTAAATAAATGGGAGAAATCACATCACGAAAATCTGTGCCTGTCGTCAGTGTGAAAA
TGCCTTTGCTGATAATTTATCCTCTAAACAAATGAGTAAATCCACAGGCAAGCAACCATATGTCTGTA
ATTGCTGTGCACTCTCATTACAGTAAGCACCAATTTTGGTGTGTGCAAGAAAATTCATTATAAGGTAAC
TGATAAAAACAGGAAATATGTGAAAATATTTTTATTAGGTGGATGAGGCCTCTTGAACAATTCAGAC
ATTCATAGTGGAGAAGTTATTCAATGAAAACCTCATGAGAAATCCTTTTCTTAATACAGCAGCACTTCTA
TAATAGATCAGAATTCACATGGTGTAGAATCTCAATGACATGAATGGAGGTAGTCTCAGTAAATTA
CTCATTCCCTTAGTCAATACCAGCATTTTTCCAGTGAGAAAATCTTGTGACAGGATAGTGGAAAACCT
TCAGGCAGCTTTTATGTCAAAAAGTGAGACAGGGATGAAAACCTTAAAAAGCCATTGATGAGATGTAT
AGCTGGGGGACAAAACATAAAGCCATCAAGCACGTGCTTGAGAAAAAATTATAATTTTGAATAAAGAC
TTTCTACTTAAATATGTGGTTGAAATGTACAATCTGAAATAACCTGGGAATATTGAATGCAGAATT
ATGTAAGAAGTAATAAGATTAATTAAGTACTGTCAAAAATA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001042415.3
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	169270
MW:	38.8