

Product datasheet for SC211638

LUC7L (NM_018032) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LUC7L (NM_018032) Human 3' UTR Clone
Symbol:	LUC7L
Synonyms:	hLuc7B1; Luc7; LUC7B1; SR+89
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018032
Insert Size:	1004 bp
Insert Sequence:	>SC211638 3'UTR clone of NM_018032 The sequence shown below is from the reference sequence of NM_018032. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TCCCGGGACCGAAGTGCAGAAATACAAGTAACTACTCTGACTCCTTCGGTAGCTGCAACCAGGAGTGAGC
CCTTCTCTGTGTTCCAGGGTCTGCTGAGGGCCGTGTCTGGTGGGGATGGGGCTGGGCTCACCTCAGG
AGTAGGGCTGGGGAGTCGTGAACGGGACTCAGGTGTGGGAAGAGGCGAGAGGGCTGTGGAGGAGCTCGC
ACGGCGCCAGGTGATGGGCTGCACAGGCACTGTCCCCTGCCTGCGTCCTGGGGCCTGTGCACTGTTGCG
TCCATGCTCAGAGTGGCTGAGACTTGTGTCTGACCAGGCCCTGCTTACCTCTGTTTTGGTTTTGTTT
TTGATATTTTTTTTCCATTGTGTTTTTACGTAGTGTCATGTTCTGTGCATATAGTGTGATTCTCCT
TTGCACTGTTTATGTTACAGTGAAGGCTCTCCTTATTAATAAATCTTCGAAAGGTCACCTTTTTAATGGC
TATCTAACTCCATATGTGGTGGGCAAGTCTGGTTGGCCTCCGGGGGGTTTCCAGGTATAGGGGATGT
AGGGCCTTGCCTGGCCTGGTCTGCGGGTCCATCGTCAGTGCCTGAGCGGCCAGCAGAAGGGGGCAGCA
GCCTCCACTGAGCCTCTGGTTCCCATTTCCAGGTTCTCCAGAGAGCGGGCATCCAGAGAGGAGTCCTG
GGAGAGCGGGCGGAGCGAGCGAGGGCCCCGACTGGAGGCTTGAGAGCTCCAACGGGAAGATGGCTTC
ACGGAGGTCAGAAGAGAAGGAGGCGGCGAGATCTGAACCCGTCCTCCGGGTGCTGTAATAGTCTGAT
AAACGTTACACAGTCTAAAATTACCCTTTATATTTTGCTGAATACAACATCATCTTTGTAGTTTAAAT
TTCTATTGTTTTGGAGTAGCTGTGAGTTTCTAGAAGGTACAGAGTTGCTCCTGTGTTCCGGGTCAT
GTTGAGTAGGAATAAATAAATCTGATGCTGCCTCCTGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_018032.5</u>
Summary:	The LUC7L gene may represent a mammalian heterochromatic gene, encoding a putative RNA-binding protein similar to the yeast Luc7p subunit of the U1 snRNP splicing complex that is normally required for 5-prime splice site selection (Tufarelli et al., 2001 [PubMed 11170747]).[supplied by OMIM, Mar 2008]
Locus ID:	55692
MW:	36.6