

Product datasheet for **SC211644**

LETMD1 (NM_001024668) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LETMD1 (NM_001024668) Human 3' UTR Clone
Symbol:	LETMD1
Synonyms:	1110019O13Rik; DKFzp586A011; HCCR-1; HCCR-2; HCCR1; HCRR-2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001024668
Insert Size:	1024 bp
Insert Sequence:	>SC211644 3'UTR clone of NM_001024668 The sequence shown below is from the reference sequence of NM_001024668. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TCCACCAACTACCTTGGGACAAGGCGCTGAATGAACCATGGAGCGGATGGCATTGTCCTGCAGTCGTAT
AGTATAGCAGTGCAGGAACAAACAGCACTTGCCAGCAAAGTCTGTGTGTACTGTTAAGTGTGTGGGAGG
CAGAGAGAGGAGCAGGGGCCATGGGCTTCACAGCATGGCACACATGTGGAACTGCAGACATTCCTCTC
ACAGCTAGAACTGAAACAAACCCTCTTGCTAGGGGTGGTCCGTGTGAGGTGTCATCCTGTCCCCCTCAT
AATTACTAATAGCTGGAAGTGGCAGCAGCCTCTACTGGGCTTTTACTGTGATGTGTTCAAGTTCATGTCC
TAGGAAGTCAGCTTTTGCCCCAGGTGGGAATCCTTATTTGGCTTAGGACTGATCCACTTCCATGTTACT
TACATCTGTGGGTTTTTGTGTTGCTGTGTAGAAAATTTGTGGCTGGTGAAAACAGCACTCCTTTGGCTG
GAGCACTTGTGTCCGTGCATGTACTTGGGTGTTTCCCTCCATCCTTTCTGATATGACCAAAAATCAAGT
TGTTTTGTTTTTGTACCTTCACTGGCATGGGCTAACCACTTCTTTTCAAACCCTCTGAACACCTTT
TTCTGATGGGTAAGTGCAGGAATATTCTATTGAAAAAGATAACAGGAAGTACAAGTGTCTTGTACCC
CTTCTCAATGTTTCTAGCCTTCACTCTCCATTGTCTTTTCTGGGCTGTATTACAGCCCTCTGTGGATC
TTCAACTCTGCTGCCTCCACTGTGATGCAGCAGTCCAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGT
CATGGATCACACCTGTAAGGTACTAATTACTGCCAGCCTGGGGAGATCAGGAGAGGTCTGCATAGTTA
GTAAGTTGGGTTTAGCTTTTGTGTGTGCATCAGTGAAGTCTGTAACTTATTGTAAATGCA
TGAAGCACTGTTTTAAACCAAGTAAAGACTGCTTGAACCTGTTGATGGAAATGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001024668.1</u>
Summary:	This gene encodes a mitochondrial outer membrane protein. It has a potential role in tumorigenesis, which may result from negative regulation of the p53 tumor suppressor gene. Alternatively spliced transcript variants have been noted for this gene. [provided by RefSeq, Aug 2011]
Locus ID:	25875
MW:	38.1