

Product datasheet for **SC211670**

NLRC5 (NM_032206) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NLRC5 (NM_032206) Human 3' UTR Clone
Symbol:	NLRC5
Synonyms:	CLR16.1; NOD4; NOD27
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032206
Insert Size:	1022 bp
Insert Sequence:	>SC211670 3'UTR clone of NM_032206 The sequence shown below is from the reference sequence of NM_032206. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AACCAGCCCCAGGCCCTTGGGGTACTTGATGGCCCCCTCAAGACCTTTGGAATCCAGCCAAGTGATGC
ACCCAAATGATCCACCTTTCGCCCACTGGGATAATTGACTCAGGAAAGAAGAGCCTCGGCAGGGCGCTC
TGCACTCCACCCAGGAGGAAGGATACGTGTGCTCCTGCTGCAGTCCTCAGGGAGAACTTTTTTGGGAACC
AGGAGCTGGGTCTGGACAAAGGAGTACCCTGCATTACGTGGGATATGTGTGATCAATTGGGGACATGCG
ACACACAATGAGGGTGTGATGACAATGCATGACACGTACGGTTATATGTGGCAGTGTGACCCCTTGACA
TGTGGCGTTACATGAAAGTCAGTGTGGCACGTGTTCTGTGGCATGGGTGCTGGCATCCCAAGTAGCAGG
ATACATGATTGTTGGTCTATATATGACACATGACAAATGTCCATGTCACAGGACTCATGGCTGGCCAGA
TGACCTCAGGCTGGCCCAAGATCTAATTTATTAATTTTTAAAGCAAATACATATTTATAGATTGTGTGT
ATGGAGCAGCTAAGTCAGGAAAAGTCTTCGCCCGAGCTGGGAGGGGAGAGTGTCCATGCACTGACCAG
TCCAGGGGCTCAAGGGCCAGGGCTCTGGAACAAGCCAGGACTCAGCCATTAAGTCCCCTCCTGCCTCA
ATCCTCAGCCTACCCATCTATAAACTTGATGACTCCTCCCTTACTTACATACTAGCTTCCAAGGACAGG
TGGAGGTAGGGCCAGCCTGGCGGGAGTGGAGAAGCCAGTCTGTCTATGTAAGGGACAAAGCCAGGTC
TAATGGTACTGGGTAGGGGGCACTGCCAAGACAATAAGCTAGGCTACTGGGTCCAGCTACTTTTGGT
GGGATTACAGTGAGTCTCCATGCACTTCACATGTTACCCAGTGTCTTGTACTTCCAAGGAGAACCAA
GAATGGCTCTGTCACTCGAAGCCAGGCTTGATCAATAAACACAATGGTATTCCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_032206.5</u>
Summary:	This gene encodes a member of the caspase recruitment domain-containing NLR family. This gene plays a role in cytokine response and antiviral immunity through its inhibition of NF-kappa-B activation and negative regulation of type I interferon signaling pathways. [provided by RefSeq, Oct 2011]
Locus ID:	84166
MW:	37.7