

Product datasheet for SC211685

DSN1 (NM_024918) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DSN1 (NM_024918) Human 3' UTR Clone
Symbol:	DSN1
Synonyms:	C20orf172; dj469A13.2; hKNL-3; KNL3; MIS13
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024918
Insert Size:	1023 bp
Insert Sequence:	>SC211685 3'UTR clone of NM_024918 The sequence shown below is from the reference sequence of NM_024918. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CATGGATCTGGATCTGGATCTTGTCTAGTCTTTATGAGAGTTTCTGCCACAAGGTGCCCAAGAGGAGA
GGAATGGGAAGAGTGCCCGAGCAGTGGTGAAGTCTGCTCGTTGCTTTGAAGATACTGG
CAGGACTGACTGTAGAACACTTTGACTTTTTTCAAAAAGTGATGGAATTTGTACATCCAAATGAATATT
GTATAGACAATTTCCAGGAATGTGAAAATGCTTGAAAGTTCAAACCTCTTTTTTGAATGATCTTC
AGATCCAGTGGCCATTCTTTATCTTTATCCTGTGAAGGTGTTTTTCAGGTTTTGAAACAATCCAAAA
ATCATTTAGGACCAAGTCTAAGGAAACATTTTAGTGGCCAAGTTGGATTCCGATTGTAAAGGAATGATA
CTAATTTTCTAGCATGGCTCTGAAGGTGATTTAGGTAGAAGAGTTTGAAGGCTGGGCGCAATGGCTCA
CGCCTGTAATCCTAGCATTTTGGGTGACTGAGGCGGGTGGATTGCTTGAGCCAGAGTTGAAGACCAG
CCTGAGAAATAAGGTGAAACCTGTCTACAAAAATACAAAAAGTTAGCTGGGTGTGGTGGCGTGTGCC
TGTAAGTCTAGCTACTCAGAAGGCTGAGGTGGGAGGATTGTTTGAAGCCAGGAGTTGAGGCTGCAGTG
AGTTCTAATTGCGCCACTGCACTCCAGCCTGAGCGACAGAGTGAGACACTGTCTTAAAAAAATTTAAAA
ATTGTAATAAATGAAAAAAGTTTGTGACATTATTTGCATCATTGGGATACATATGTCACCTTCACA
AGATGTTCAATTTGAAGGAAATACCACTCATTCTCTATGCTGTTGTCTGTAGTGTGCTTCAGTTTTT
CATATTGAGTTGACCTAAATCCTGGATTATGACAAGAAAGGAGTAAGTACTACTATTATTGTTCTAT
TTGTTTATAATCTGTATTATAAATTGCACATAATTTAAAGCTTTCCCTTGTCTTCA
ACGCGTAAAGCGCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCGCCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_024918.4</u>
Summary:	This gene encodes a kinetochore protein that functions as part of the minichromosome instability-12 centromere complex. The encoded protein is required for proper kinetochore assembly and progression through the cell cycle. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2009]
Locus ID:	79980
MW:	38.8