

Product datasheet for SC211712

EHBP1 (NM_001142615) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	EHBP1 (NM_001142615) Human 3' UTR Clone
Symbol:	EHBP1
Synonyms:	HPC12; NACSIN
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001142615
Insert Size:	1017 bp
Insert Sequence:	>SC211712 3'UTR clone of NM_001142615 The sequence shown below is from the reference sequence of NM_001142615. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AAGAAAGAGGAGAAATGTGTTCTTCAGTAGCCATCAGATCAGAAAGAATCTCTCCCAACATTTTAGAGT
CTTGCTTCCCAACCAGAAAAAGTCAGACTCATTGTTGATTTAAACTTTTAACATTTTGTGGCTGG
ATTGTACTACTTTACCTCTACTTTACCACCACCACCTTTTCTCCCTCCTTTCCAAATAATATACAGA
ACTCCAAAATAGCTTCATTTAAGGATTTTTTGTGAGTTAACAATTTCTTGAAATCCTGTGAAATAGA
TTTGCACAGACACCTTGTGAGTGATTGGTATTGGAGGTGTCAAGAAACTGTTGAAAAAGAACAAAAA
CACTTCCCTCGTTATTTTCTCTCATTTTTTGTGAGAGGAAAAATTGAAACATTATTCTTGTGTTGTT
GGTAATAGCATAATGACAGTGGGAGGGGTACAAGGGGATAAGAAAAATGTCATGATTTTTTCCGGTCC
TGCCACATGTAACACTTACTCTGTTACCTAAATTTTATAGTTAGATCATATCCAATCTACTTATTAAC
TGTGTTCTATTTACCAGTGGAGTTTTTCTGCAGTGGTTGCGTTTCACTGTAAGGATAATGGAGTTCCTC
TCCTCTGCTTTCTCAGAGGATGGTCTTTAACATAGCCAGAAACAAGCCCTGTGGTTTGAAGGTGAGC
TGTGAGGATGGGACTAATTGATATGCACCAGTTTACAAAGACAGTCTTATCATCCGAGAATACACCATC
TTTTTCTCTGGATAATTATTTCTTACATCATGCTTGATTCTACATTTTGTGGTCTCAACATTGGCT
CACGAATGCTGTTAATATTTTCTGTATTGATAAAAAGTCTGTCTTGCCACTACAAGTAAATCCCCCA
TTTAATATTTTCTTCTTAGCATAGCACTGTCATTTTTGTGAAATGGTTATGTTTATTATTACAAT
ACTGAGTCATATATAATTTTCAATAAAAGCAGAACTTTCTTACCTTAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001142615.3</u>
Summary:	This gene encodes an Eps15 homology domain binding protein. The encoded protein may play a role in endocytic trafficking. A single nucleotide polymorphism in this gene is associated with an aggressive form of prostate cancer. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2010]
Locus ID:	23301
MW:	40