

Product datasheet for SC211719

CCDC62 (NM_201435) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CCDC62 (NM_201435) Human 3' UTR Clone
Symbol:	CCDC62
Synonyms:	CT109; ERAP75; TSP-NY
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_201435
Insert Size:	1177 bp
Insert Sequence:	>SC211719 3'UTR clone of NM_201435 The sequence shown below is from the reference sequence of NM_201435. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GCTCAAAAAAATACCTTTGTCAGTTATAGAGGAAACAAAAGGCAACTTCAGTATTCATCGTGATCAGC
AATTTCTCATCTATGTGGAAGGCAGAAAGCAGACCAATACTGAATGAATACTTAACCGTAAACTGA
AAGAGGATTCTAGTTCTTCATAAACGGCACTTAATTCCAGCTGGGAGCAGAACTAGAAAGTTAATTTTT
AAACATCTACACTTCATTTTCAAGTTAACCATTTTTGTGCTGAAGAAATATTTTCATGTGAAGAAAGT
AGACCTTATTGTACATATAGAAAGTTGGAATTATGCTAAGAATGAAAAGACTTCTCTGTAAAGATACG
GACTACAGTTAAATGCTAGAGAAGCTCTTTAAAAATGTGAATGTCAAATAGAGAAAGAACCCTGCATAG
AAAGTGCTGTTTTAACTATCTGATTTTTAAAAATCTGTGCATACATTTAAATTCTAAACAATAGCTTA
TCAGAGTCAGCTCAAAATATATGAGAAACAGTATTCTCTCATGGTTTTAGCTTTTGACTTTGCTGTGTA
AATAGACATAAGGTGCTTTGATATAAAATATAAAATGTAAGTGGAAATAGCTCGAGGTCCTTCTGTCC
CAAGCTGAGCAGAGCCCCATCTTTCTGGGTCTATATTAGTCCACCTACTGACACAAACAAAAGCTTGC
TGGAAGATCGAGTTTTAGACGCATTTTTAAAAATCTTAAAGACTAAAACACTTCCATTTTAACTTGTA
AGTAATTTAATTTTTAAAGATTATACTATATGCCTCTGTGCTTCTCTAAAAGAATAGATCAACTTCA
GTCCATAAAAGATATTTTTAATATTAAAGAAAAATATGTTTCCTTGTTTCTTTTATTTTACAGGAG
TAAATAAGGAAGGAACGTTTCATCACTTTAACTGAACCTGGCAAGTTAATTTCTCAGGAATGGGGAT
GTATTTTTTTAAGCATTGCAGATATCAAAGTTCTATTGTGCTGAATAAATGCCCTTTGTTAACAGGTC
TCAGTGTTACTTCAGAGATGATTTTCATGAGCGGGAAGAGCCACCTCATTTAGCCTTTTAAATAGGG
TGCTGAAAAGAGAAGTTGGGTGAGGTTTGGTTATGTTCTTTATCAATGCGTGCTTTCTCCCTCAGGCTG
GGCT
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_201435.5
Summary:	Nuclear receptor coactivator that can enhance preferentially estrogen receptors ESR1 and ESR2 transactivation. Modulates also progesterone/PGR, glucocorticoid/NR3C1 and androgen/AR receptors transactivation, although at lower level; little effect on vitamin D receptor/VDR.[UniProtKB/Swiss-Prot Function]
Locus ID:	84660
MW:	46.5