

Product datasheet for SC211731

SWT1 (NM_017673) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SWT1 (NM_017673) Human 3' UTR Clone
Symbol:	SWT1
Synonyms:	C1orf26; HsSwt1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_017673
Insert Size:	1008 bp
Insert Sequence:	>SC211731 3'UTR clone of NM_017673 The sequence shown below is from the reference sequence of NM_017673. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TGTGAAGACATGCTCAACTATAGGATAAGTACTGATTGTAACTTTAAAGGAATTGCATTGTCTT
AAGAATAACAGAGTAGTTTTCAATCTGGTCACTCTTTGGGCCAAACCAAGAGAATTTAAGAAATGT
TTCATAGGTATAAAAAGGTGATCGCCTATTACTGACAGTCTCATTGTAGCTCTAAAAAGCCTAATGTAT
CCACTGTGGAATAAACTCCATAGACTCAACTCTTCTGGAGTTCACAATGCCTCCCTACCTCTATTCTTG
ATAGTGAGCCAGTATCTCTAAGGAAAGGAGAGCCGTATTTTAGTCATCCTAGGGCAGGAGCCACTGTG
GCATCTGAGGAAGACCCAGGGTACAGTATTAGCATTTCTTTCAACCCTAAGTATGGGCATGTGACCAGA
GACACTGCCAAACCCTAAGAAATCATTAGGCATGATCATGCTCACCATATTTCACAAAATGAAATCTTA
ATGTGGGCTGTTTTTATGATCAAACTGTCATCCATGTTGACCTGATGAAAATAACACCCTCCTTCACC
TCCCTCAGGGTAATGAAAGTAGAAATTTGACTGAAGCTCCTACTTCTTATACAGCCTTCTCTGTG
GGCTTGTCACCTTTTCAGTACATTTCTAGTGGGGAAGAGCAGGGCGACCCCTCTGCTGGGTTTCTCCTTA
AATGTGAAGCAAAATAAATTTATATAATGGAATTTGGTGGCAGTACTGATGTTTGTCTAGCTAAAA
CCAAGTGTAGTTGACTTTTATAGTATAGGATAAGTTGGAAGTCCTTGTGTTACTGGAATATATTATGC
TTTTGTTAACACATTAGATTAAATAGAATTTGGTATTCTATTATAAATGTTTATTTTTTAAAGCAATAA
ATAACTTATTTTTGGTAAGATTTGGCTTTTACTTCAAATTGTTTCTGATATGTACATAAGGTATGAT
TATAAAAGTGACTGTCTGCAATAAAAGAGAAGTTGTTTCAGGA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_017673.7</u>
Locus ID:	54823
MW:	38.8