

Product datasheet for SC211735

BORA (NM_024808) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	BORA (NM_024808) Human 3' UTR Clone
Symbol:	BORA
Synonyms:	C13orf34
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024808
Insert Size:	1010 bp
Insert Sequence:	>SC211735 3'UTR clone of NM_024808 The sequence shown below is from the reference sequence of NM_024808. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GCAAGCCCTTTTCAATGCAGCAGTCCAATGCCTCTGTGAGAACTAAAGCTTAAGAGTTT
CTCGCATATATCGTTGTGCACAGGATCAACATGATGGTGACTGGGAAAAAATTACTTCAAGTAACATGC
TTAGCTTTCCCTCCTTAATGTGAAAAATCAAGGGCTTACTGACATAGGAACAACAGAAATGCTCCTGGA
ACTTCAAGTTGCTGAATTATAAGTTTATTTTTATCAATAAATATTTTATACTTACATTGAGTGATGT
GTTTAACAACAAATTGTGACAGAGCTGAGTGCTCCTATCTTACAGGGTCAATGAACTACTTATTAAGCC
TTACTGGTAGCACTGAATTTAGCAGTTCTGAGAACATGTGAAACTATGTTAAACTGAAGGCACTATAT
ATTTTACATAAAAGCTTGAACATACAGATGAATTATAACCTATGTGAAGAAATCTTAGATATAAACT
AACTTTTCAAAGATACAAAAGAAATTAACAGGTTTCTGAAATTTAGTTCTTGGTCTGTTACCTCT
GTGGGGAAAAATCTTAGTTCCAGTGATAACTGTTCTAGTTACTACTTTTAAAGTATGTAATACTAGAAA
GGTAGTACTAGTGACATCATCACGTGTATTGTTATCTATGGGGCAATGTGTGGTGCCGAGAATAAAAT
ATACCTCATGCCTAGGTAGGGACATCCTTTCCAGCTCAAACGTGGGTAGGGATGTGGGAGAATAAGAA
TGTGGGAGAACCAAGAGAAAAAGTGGGGCTGGGAGAGTGGAGTTCCCGTAGGGCATAGGCCTGTGAAGT
AACACTGGGGCAGATATGTATGTTATATACAACATTTTTTTTAAAAAACTTATATCCATGTTGGGAGTA
GATGGGTATATAACAGTTTGAAATACTATCTTTGGAGAATGTATTTTGTATTATATAAATCAACTTTT
AAAACTGTCTCATTCAAAGGGAATAAAGACCTGTGTTATCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_024808.5</u>
Summary:	BORA is an activator of the protein kinase Aurora A (AURKA; MIM 603072), which is required for centrosome maturation, spindle assembly, and asymmetric protein localization during mitosis (Hutterer et al., 2006 [PubMed 16890155]).[supplied by OMIM, Mar 2008]
Locus ID:	79866
MW:	38.2