

Product datasheet for SC211738

CEP128 (NM_152446) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CEP128 (NM_152446) Human 3' UTR Clone
Symbol:	CEP128
Synonyms:	C14orf61; C14orf145; LEDP/132
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152446
Insert Size:	1030 bp
Insert Sequence:	>SC211738 3'UTR clone of NM_152446 The sequence shown below is from the reference sequence of NM_152446. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CAACCCAAAAAGAGGAATATGGGAGCTAAAAAGCAAATGTAATTTGTTATTTTACATGAGTATGTTA
CAAATAATAACATCTCTATTCTTACAGCAATTTGGCCAGATTATCTAACAGACATACCTGCAGCTTTG
GCTCTTTGGTATTGCCAAACATTGACAAAAGTGACAATACTGTTGGTCCTTGTGAATGGTAAACCAATC
CAAATAATATCAGATCATGAATGATGTGCAGCTAATTTATTTGCAACTTATACTAGCACAGAAAGAATG
ATATTTTTTCACTTTGTATTTTCTTTTCAAGTGTGGTCAGTCACATTTTCCCAACCAAAATTTGGAT
TTCTAATTGAGGCCTTTTAAAAATTACTGGCTAAAAAGAAGGAAAATGAGTTTCCCTTGAGGCAAAAGT
TGCCAGCACTTTAGATTCAGGGAGCATTATATACAACCTGATTTTTTGAAAAATTTCTCCCTAATAAA
TTATCATGTGTTATAACATGTGGAGATGTGCATATATAACATTAGACAGTCCACTTTATAAATGCTCAG
GACTTTTATATCTATCATCTTAACAGTCCACAACCTGCATACCGTAGCTGAAAATTTTCCAATCATCC
ATTTTAAATAGCAATAAAGAATTCATAAATGATATCTGAACTGTAGAGCACTGTAAGAATTCATAGCA
CAATTTATTTATTTGTCTTTTGTGGCCAAATATCCTACTGTGATCACTTCAGTCATATCCTGGACCGT
CTTTGGTTTACTGTGACTAAAATCATAGCACCAAACTCCTTTGTCCATGAGATTTATATAATCTTGT
ATAATTTGGTCCTATCTGCTTCTACCATCCCAACTTTTTTCATCAAAGTGTCTAAGAACGATTCCCTT
TACACGCTGATCTCCTAGTGGAGAAGGATTTATATTAATGCTAAGTAATAGTAACTTTTCAAGAACTG
AATCCTATTGGTCATGGATTCTGCATAGACAAGGGCAATAAAGCAATACCTTTGTAATATGCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_152446.5</u>
Locus ID:	145508
MW:	40.5