

Product datasheet for SC211745

MRFAP1 (NM_033296) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MRFAP1 (NM_033296) Human 3' UTR Clone
Symbol:	MRFAP1
Synonyms:	PAM14; PGR1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_033296
Insert Size:	1042 bp
Insert Sequence:	>SC211745 3'UTR clone of NM_033296 The sequence shown below is from the reference sequence of NM_033296. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CGGCGGATAGAGAAGAGCGAGTCGTCGTACGCGCGGTTCAGCCAATGGATTCTGGTCAAC
TGGTGGAGATTGGCTGACACCTGGAGAAGCCGAAACCAGAGAGCCTTTGTTTTCTCTTTTCTGT
CTATGCTCTGTCTCACTTAACACTACGTTTTCTGCTATGGTCTGTGGTTGATGACCTCAATATGAGTTT
CGATTGTTAACGTGTTTTGTTGGGAAGTAATTTGTTTGAAAATGCTCTCACATACAGGAATTAGGG
CCTAGATTGTAAGCTCTTGACAGCAGTCACATTTGTTCCCGGGCTTTGGTGGTTATTTCTAAATTTTGA
GGTGCTTTGCTATTTCTGTGTGACCTGATAGCTCCCTGGAACCTTTGGGTCTGTGTGTGACACATGAGA
CTCACAGTTGGAGTTCTCCAGCTCTGGAGGTGCTGAAGGAGCTGCATTAATTCTGGAAGACGACTCCAT
GCAGCAACTACTGAAGAAAGGACCAGACTTCAACGGGGAGTGTGGATGGGCCGACCTGGCTGGGACTCG
TGAATCTGGAGAAGAGCTGGAGAATGGATAGTATTGCTGTATTGAGACTTTAATTTCTGTGTGAGA
CCAAAGGAGGAGAGATGTGTTTTGTTCAAAATTTAAATTTGTTGTGGTACACTATCTTATGTAACCTGT
CTGGTGAGTTTGTGTTGGACAACCTAACTCAGCTTTATTTGACATGGAACCTAAAAAGAGATAAGATC
TTGATATTCTGTACAAGTTGATGTAATACCTGATGCGTTTTAGAGGACTTGGCATAAAATGAAAGATT
GGCAAAGGCCCTTGAGGGCTTGGGGATGACAGTATGGAACCTGTCTGCATTGGACCCTAAACTGGACTA
GAAGAGGCATCTTCAAGTTTCATACGTTGTCCAGCTGTAAGTTTCATTTGAGTAGCAGACCTAACAAATA
TTTGAGGTCAGAACCTACCATGTTAAACAAACAAAACCTTACCATGTTAATAAAGTATTTCATTTGC
TTGAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_033296.3</u>
Summary:	This gene encodes an intracellular protein that interacts with members of the MORF4/MRG (mortality factor on chromosome 4/MORF4 related gene) family and the tumor suppressor Rb (retinoblastoma protein.) The protein may play a role in senescence, cell growth and immortalization. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013]
Locus ID:	93621
MW:	39.4