

Product datasheet for SC211748

CRMP1 (NM_001313) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CRMP1 (NM_001313) Human 3' UTR Clone
Symbol:	CRMP1
Synonyms:	CRMP-1; DPYSL1; DRP-1; DRP1; ULIP-3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001313
Insert Size:	1026 bp
Insert Sequence:	>SC211748 3'UTR clone of NM_001313 The sequence shown below is from the reference sequence of NM_001313. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GGCCGCTCCAACATCACCAGCCTCGGTACGTGGATGCGCGGAGGAGCTAGCCTGAAGGATTCTGGG
AATCATGTCCATCCCTTTTCTGTCAGTGTGTTTGAACCCACAGTTTTAGTTGGTGTGATGGAGGGA
GGGGGAAGTCGAAGGATGCTCTTCCCTTTTCTGTTTAGGAAGAAGTGGTACTAGTGTGGTGTGTTTGC
TTGAAATTCTTGCCCCACAGTTGTGTTTCATGCTGAATCCACCTCGGAGCATGGTGTGTTTTCATTCCCC
CTTCCTAGTGAACCACAGTTTTAGCATTGTCTTGTCTGTCCTTCCACTTCTAACTCCACTGGCTCC
ATGATTCTCTGAGTGGTGGTTCCTTGCACCTGTAGATGTTCTAGGATAGTTGATGCATGTTACTAAA
TTACGTATGCAAGTCTGTGAGTGCCTGTGAGGGGACATCGCCAAGGACTGACTGAGACACGATGCCGAG
ACCTCAAGCCCTGAGGGGAGTCCCAAAACCTTACAGTGAAGATGTTTACTCATTGCCCCACCTCTG
GTCCACACTAGAAAGAAGCTCGCCCCACCTCCACCTGTGAGATCCGTGAATTCTCGGAATGGCAGGGGA
AGCCTTGCACTAGGTTGCAGAGAAGCATCCTCCACATCCTGTGTCAGAAACCTGGTCTCCGTGGCACT
TGTAACCTACCGTGCTGTCTTCTGGTCTGTGTGTTCTTCAAGCCAGCTCTAGGCTTCAGGCCGAGCC
AGGTTCACTCAGAAAGAGGTCTCCCCATCCCCATTGCGGGCTGACGATGGGGGCTGATGGCTGCCC
CTGCGTGGCTGAGTCTGGTCCCTCTGAGGCAGTTGACGGGGCAGTCAGATTTTAAAGTTTTGTACA
AAGTTTTCTTTGTAATCACTCCATTTTTACTTAACAACCAACTGTTGTGGCTCTTATTCTGAATT
CAAAGCTTGTGAAAAATAAAGAAATGAAGTCCCACTGACCTGACTAGGGCACTTTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001313.5
Summary:	This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008]
Locus ID:	1400
MW:	37.2