

Product datasheet for SC211749

CRMP1 (NM_001014809) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CRMP1 (NM_001014809) Human 3' UTR Clone
Symbol:	CRMP1
Synonyms:	CRMP-1; DPYSL1; DRP-1; DRP1; ULIP-3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001014809
Insert Size:	1026 bp
Insert Sequence:	>SC211749 3'UTR clone of NM_001014809 The sequence shown below is from the reference sequence of NM_001014809. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GGCCGCTCCAACATCACCAGCCTCGGTACGTGGATGCGCGGAGGAGCTAGCCTGAAGGATTCTGGG
AATCATGTCCATCCCTTTTCTGTCAAGTGTGTTTGAACCCACAGTTTGTGTTGGTGTGATGGAGGGA
GGGGGAAGTCGAAGGATGCTCTTCCCTTTTCTGTTTAGGAAGAAGTGGTACTAGTGTGGTGTGTTTGC
TTGAAATTCCTTGCCCCACAGTTGTGTTTCAATCCACCTCGGAGCATGGTGTGTTTCAATCCCTCC
CTTCCTAGTGAACACAGGTTTTAGCATTGTCTTGTCTGTCCTTCCACTTCTAACTCCACTGGCTCC
ATGATTCTCTGAGTGGTGGTTCCTTGCACCTGTAGATGTTCTAGGATAGTTGATGCATGTTACTAAA
TTACGTATGCAAGTCTGTGAGTGCCTGTGAGGGGACATCGCCAAGGACTGACTGAGACACGATGCCGAG
ACCTCAAGCCCTGAGGGGAGTCCCAAAACCTTACAGTGAAGATGTTTACTCATTGCCCCACCTCTG
GTCCACACTAGAAAGAAGCTCGCCCCACCTCCACCTGTGAGATCCGTGAATTCGGAATGGCAGGGGA
AGCCTTGCACTAGGTTGCAGAGAAGCATCCTCCACATCCTGTGTCAGAAACCTGGTCTCCGTGGCACT
TGTAACCTACCGTGCTGTCTTCTGGTCTGTGTGTTCTTCAAGCCAGCTCTAGGCTTCAGGCCGAGCC
AGGTTCACTCAGAAAGAGGTCTCCCCATCCCATTCGGGGCTGACGATGGGGGCTGATGGCTGCC
CTGCGTGGCTGAGTCTGGTCCCTCTGAGGCAGTTGACGGGGCAGTCAGATTTTAAAGTTTGTACA
AAGTTTCTCTTGTAACTCACTCCATTTTACTTAACAACCAACTGTTGTGGCTCTTATTCTGAATT
CAAAGCTTGTGAAAAATAAAGAAATGAAGTCCCACTGACCTGACTAGGGCACTTTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001014809.3</u>
Summary:	This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008]
Locus ID:	1400
MW:	37.2