

## Product datasheet for **SC211753**

### SCMH1 (NM\_001031694) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** SCMH1 (NM\_001031694) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** SCMH1  
**Synonyms:** Scml3  
**ACCN:** NM\_001031694  
**Insert Size:** 1022 bp  
**Insert Sequence:** >SC211753 3'UTR clone of NM\_001031694  
 The sequence shown below is from the reference sequence of NM\_001031694. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATTGACCGGCTGAAGCAGGGCAAGTTCGACCAGGAGAGCAGCCTAGACAACCAAGTGGCAGCAGGT
GGGGGCATTCTTCTAGGAATGAGGGGCATCAGCCACCCAGGCACCTCAGTGGGTTCCGGGCCACCT
CAGGACTCCAAGAGGCTGTGTGGAGCCACCACTCCTAGCCACAGCTGCCATGATAAGTCCTTCCATGAA
GGACTGAGGAGGAGAGTGGGGTCCAGGGCTGGTGTCTTCCCTCAGCTCTGCCGGGGCTCTAAG
GTCCCTCTATTTATTTCTCAACCCTGGCTGGCCTCTCACCAGGAGTTTAGGCTGAATGCCTTCCACGTG
ATGGAGGAAAAGGCCAACTCTGTCTGTCTTGTGTGGCACCCCATCGCCCCACAGCTCGTACCTTCT
CACCAGATTCCTGAATCCAACTCGTGGTGCAAACCTCTACCTTTTTTACAAAAGATCTTATTGTT
AATTTATTGTTTCTGGCACTTGGCAAACCTGTAGTTAATACTCTCCACACTAGACACTGGGTTTC
AGGAGGAGGGAGACTGCCCTGCTTTGGTCCCAGAGAGGCCCTCTGCAGATAGGCGTGGCCCTCTCAG
AGGACACTACCCTAGGGCACTTTCTCTTTGAGGTGGAGAGACCCATAAAGCCTTGACCACATCACTCCA
TATGGGGAGGAGAAGGATCCCTGTACCTTCTCCTCTCTTACGGGGCCCTTTTGCAGCCCTAGGCCTC
ATCTGTGGGAAGGGAGTCCCTGGCTTATACTGCCCCACCAAGCTCCTTGCCTGGCCAGAAGTCTG
TCGAAGAAAATCAGGCCGGAAGGCCAAGAAGGCCTAAGGGGGATGGGAGGGCAGGTTTTCCAGGCTGG
AGTCGGTTCACCCACTCGCCTGTCCACAGGCTCCTTGTAAAGCAAGTCAGCAGCACAGCTACTCACGC
TGCCATCTGGACTATTTTATGTCAATCTGTTTATAAATAAAAACCAATATAGATA
ACGGCTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_001031694.3</a>
<b>Summary:</b>	Associates with Polycomb group (PcG) multiprotein complexes; the complex class is required to maintain the transcriptionally repressive state of some genes.[UniProtKB/Swiss-Prot Function]
<b>Locus ID:</b>	22955
<b>MW:</b>	37.8