

Product datasheet for **SC211757**

SLC6A12 (NM_001122848) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC6A12 (NM_001122848) Human 3' UTR Clone
Symbol:	SLC6A12
Synonyms:	BGT-1; BGT1; GAT2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001122848
Insert Size:	1017 bp
Insert Sequence:	>SC211757 3'UTR clone of NM_001122848 The sequence shown below is from the reference sequence of NM_001122848. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATAGCCGGGGAGAAGGAGACCCATTGTAGGGTGTGGCCAGAGGCCAGGCGGCTCCTAAGCCGGGAACC
TAGGTCAGGGCCACCCTCCATTCTCAGCGGACAGCCTCTGCCTCTGTCTCCTGCCACAATCCTGCTGGG
AACCTCTGGAGAGCCACAGGCACCCCACTGGAGGCCAGACTCCTCTCTTGTGCTAGCTGGAGCAGCT
CCTTCCCCTTTGCTGATAACACCACCACTGGGACGTGCCATGTTGGGACGCCACTCCCTGTGGGAAGGC
ACCATCGTTTTTATAAAGGGGGTCTTTTTGGAGGCCGCGCATCTGATTGCAACACCTCGAGTTATGAGG
ATTCCACTGTGGGGATGCCTCTTGTTAGAGCGTACTGCATTTGTACACGGGGAGAGGAGCTATAATTGG
AACGCACACTGCCGTCCAATGTGGAGAGCCTGATGGGACAATACCCTGTTGGAAGTGACAACCTGAACAC
ACTGTGTTGGATCGGAGGTTCCGTTAGGGGATCCTTCCTTAGGCTTAACGACAGAGGCAAGCCTTTGCA
TGCCGTCAGTCTGGAGTTTCTCCGAGTCTCTCATGGCATCTCCAGCTCCTGCCCTAGTTCCGCACTGT
TCTTGCACTGTTTCATCAACTCCTGGAGCATTGGAATGGAAGGGGCTTGGGAGATGATTCCTAGACTTC
ACAAACACTCGGCATGCCTCCCTGCACTGTCCGTTCTCTGCCCCAAGGCCGATATTGCTAACTGATCAC
AGATTCTTTCCACCTCACAATCCTTCCGAATGTGCTCCAGGCAGCACCATTGCCATCCTGCTTCTAA
CGCAAACCCCTGACTTTCATGGATGAGGAACCTGGAGACCAAGAGACAAAGGGGACTTTTTCAAGTTCAC
ATGGGGACCCCTTCTTGGGGGCCAGAGATATGACTAAAACCTTATCTCCTTGTGCTCAGGCCAGTGTC
TTCCCATTAACCCCTGCCTTAGTTAACAAGTGTGTATGGATTGCCTTCCA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001122848.3</u>
Summary:	Transports betaine and GABA. May have a role in regulation of GABAergic transmission in the brain through the reuptake of GABA into presynaptic terminals, as well as in osmotic regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	6539
MW:	37.5