

Product datasheet for SC211766

SEPTIN6 (NM_145802) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SEPTIN6 (NM_145802) Human 3' UTR Clone
Symbol:	SEPTIN6
Synonyms:	SEP2; SEPT2; SEPT6
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_145802
Insert Size:	1039 bp
Insert Sequence:	>SC211766 3'UTR clone of NM_145802 The sequence shown below is from the reference sequence of NM_145802. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGAGACAAAGAGAAGAAAACTTTTTTAATCCTTGCTTCCAGCAGCTGCACTAAGTCTAAAGGAGAAGA
CTGCCATTATAGAAGAGTTAGGGTTCCATATTGTCTCAAATCAGAAATCAACCAATTTCTCCTCCCTC
AAACTGCAAGCACACACATACACCACACCACTCAACAAGTGTTTCATGTGTCCCTGTGTCCAGGCAAG
AAGCTCTCTTCTGACTCACATGGTATTTAAATGGAAGTGTCTTGTCCCTAACTAACAAGGCAGGAAAA
GAACCATCAGAGCTGGAAAATGGACGAAATGTAACCTCAGAGAAACAACCTACAGGACCACTCACCCAAG
TGTAAGTGACTGGGGCAGGACACCTCAGCTGTGGGTATGAAAGTACTGTTCTGTTTACAAGGTTTTGTT
TGAGTTTTATGTTTTCTTTTAAACATTTCTCTGGTTCGATGGGTTGACTGTCTACAGCCACTGTTAAA
CATTTCTGAATATGCAAGAGAAAGTCAAGTGACATTTGTATCTTCTTCCAGCATTCCGAGACCTTCTATA
GATTCAGCAAAGGGGAAAAATGTATCCACTATCTAACACTTAGGTAGAGAAGGGAGGGGTTTAAGCT
TAGTGAGGGCAAATTAATCCATTCCACCTCCGAGACCAAGTTAGGGTTTGGAGAGAGGTTTCTGCATC
ACCTGGGATCTGGAGGGAGAGCTTTGATGTTGGTAAATCTGCCTTGAATTCATTGGTTTAACTTGCATC
AAAATACCATGTGAGTGTCTCATTCTCATATATCCCCTACCACCATCACCACCTCTGTGTTTCTGTT
GTCTCTTGAGAGAGCCTCTTTGCATGTTTTCCAGAATCTGTGTGTTTTTCTTTCTCTCTTTGTT
CTTTTTGCTCAAAGGTGTGACCAGTCATTGCCCTCTGGGGCTTTCATTCTCCAGGAGAAACATCCAG
AACCAGCACTGTTTAGCCTGATACCTTTCTAATGTCCATGTCAATTTTCAATAAAATTCAAAGAAATGC
TAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-Mlul



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_145802.4
Summary:	This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis. One version of pediatric acute myeloid leukemia is the result of a reciprocal translocation between chromosomes 11 and X, with the breakpoint associated with the genes encoding the mixed-lineage leukemia and septin 2 proteins. This gene encodes four transcript variants encoding three distinct isoforms. An additional transcript variant has been identified, but its biological validity has not been determined. [provided by RefSeq, Jul 2008]
Locus ID:	23157
MW:	38.9