

Product datasheet for **SC211799**

PCYOX1L (NM_024028) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PCYOX1L (NM_024028) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PCYOX1L
ACCN: NM_024028
Insert Size: 1029 bp
Insert Sequence: >SC211799 3'UTR clone of NM_024028

The sequence shown below is from the reference sequence of NM_024028. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTGATGCACAAGGTCAAGACTGAACTGTAGGGCTCTAGGGAGAGCCTGGGAACCTTTCATCCCCACTG
AAGATGGATCATCCACAGCAGCCAGGACTGAATAAGCCATGCTCGCCACCAGGCTTCTTTCTGACC
CCTCATGTATCAAGCATCTCCAGGTGACCTACTGTCTGCCTATATTAAGGGTCCACACGGCGGCTGCTG
CTTTTTTTTAAAGGGGAAAGTAAGAAAAGAGAAGGAAATCCAAGCCAGTATATTTGTTTTATTTATTTT
TTTTAAGAAGAAAAAAGTTCATCTTCACAAGGTGCTTCAGACTTGGTTTCTTAGCTAGAAACCAGAAGA
CTACGGGAGGGAATATAAGGCAGAGAATATGAGTCTTATTTTACTGTTTTTCACTACCTACTCCC
ACAATGGACAATCAATTGAGGCAACCTACAAGAAAACATTTACAACCAGATGGTTACAAATAAAGTAGA
AGGGAAGATCAGAAAACCTAAGAAATGATCATAGCTCCTGGTTACTGTGGACTTGATGGATTTGAAGTA
CCTAGTTCAGAACTCCCTAGTCACCATCTCCAAGCCTGTCAACATCACTGCATATTGGAGGAGATGACT
GTGGTAGGACCCAAGGAAGAGATGTGTGCCTGAATAGTCGTACCATATCTCCAAGCTTCTGGCAACC
AGTGGGAAAAGAACATGCGAGGCTGTAGGAAGAGGGAAGCTCTTCTTGGCACCTAGAGGAATTAGCC
ATTCTCTCCTTATGCAAAGATTGAGGAATGCAACAATATAAAGAAGAGAAGTCCCAGATGGTAGAGA
GCAGTCATATCTTACCCCTAGATGTTTCATCCCAGCAGAAGAAAAGAAAGGTGTTGGGGTAGGATTCTT
CAGAGGTTAGCCTGGTACTTTCTCATCAGACACTAGCTTGAAGTAAGAGGAGAATTATGCTTTTCTTTG
CTTTTTCTACAAACCCTTAAAAATCACTTGTTTTAAAAAGAAAGTAAAAGCCCTTTTCATTCC
ACGCGTAAAGCGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

| | |
|--------------------|--|
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | NM_024028.4 |
| Summary: | Probable oxidoreductase.[UniProtKB/Swiss-Prot Function] |
| Locus ID: | 78991 |
| MW: | 39.8 |