

Product datasheet for **SC211803**

P4HA1 (NM_001142596) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: P4HA1 (NM_001142596) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: P4HA1
Synonyms: P4HA
ACCN: NM_001142596
Insert Size: 1035 bp
Insert Sequence: >SC211803 3'UTR clone of NM_001142596
The sequence shown below is from the reference sequence of NM_001142596. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGACCTTGTACGTTGTCAGAATTGGAATGACAAACAGGCTTCCCTTTTTCTCCTATTGTTGTA CTCTTA
TGTGCTGATATACACATTTCTAGTCTTAACCTTCAGGAGTTTACAATTGACTAACACTCCATGATTG
ATTCAGTCATGAACCTCATCCCATGTTTCATCTGTGGACAATTGCTTACTTTGTGGGTTCTTTTAAAG
TAACACGAAATCATCATATTGCATAAAACCTTAAAGTTCTGTTGGTATCACAGAAGACAAGGCAGAGTT
TAAAGTGAGGAATTTTATATTTAAAGA ACTTTTTGGTTGGATAAAAAACATAATTTGAGCATCCAGTTTT
AGTATTTCACTACATCTCAGTTGGTGGGTGTTAAGCTAGAATGGGCTGTGTGATAGGAAACAAATGCCT
TACAGATGTGCCTAGGTGTTCTGTTTACCTAGTGTCTTACTCTGTTTTCTGGATCTGAAGACTAGTAAT
AACTAGGACACTAACTGGGTTCCATGTGATTGCCCTTTCATATGATCTTCTAAGTTGATTTTTTCTCT
CCCAAGTCTTTTTTAAAGAAAGTATACTGTATTTACCAACCCCTCTCTTTTCTTTAGCTCCTCTGT
GGTGAATTAACGTA CTGAGTTAAAATATTTGATTTTTTTTTTTTTTTTTTAAATGGAAAGTCCTGCATA
ACAACACTGGGCCTTCTTA ACTAAAATGCTCACC ACTTAGCCTGTTTTTTATCCCTTTTTTAAATGA
CAGATGATTTTGTTCAGGAATTTTGCTGTTTTTCTTAGTGCTAATACCTTGCCTCTTATTCTGCTACA
GCAGGGTGGTAATATTGGCATTCTGATTAATACTGTGCCTTAGGAGACTGGAAGTTAAAAATGTACA
AGTCTTTTCAGTGATGAGGGAATTGATTTTTTTTAAAAGTCTTTTTCTTAGAAAGCCAAAATGTTTGT
TTTTTAAGATTCTGAAATGTGTTGTGACAACAATGACCTATTTATGATCTTAAATCTTTTTTAAAAAAA
ACGCGT AAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001142596.2
Summary:	This gene encodes a component of prolyl 4-hydroxylase, a key enzyme in collagen synthesis composed of two identical alpha subunits and two beta subunits. The encoded protein is one of several different types of alpha subunits and provides the major part of the catalytic site of the active enzyme. In collagen and related proteins, prolyl 4-hydroxylase catalyzes the formation of 4-hydroxyproline that is essential to the proper three-dimensional folding of newly synthesized procollagen chains. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]
Locus ID:	5033
MW:	40.3