

Product datasheet for **SC211834**

SUMF1 (NM_001164675) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: SUMF1 (NM_001164675) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: SUMF1
Synonyms: AAPA3037; FGE; UNQ3037
ACCN: NM_001164675
Insert Size: 1028 bp
Insert Sequence: >SC211834 3'UTR clone of NM_001164675
 The sequence shown below is from the reference sequence of NM_001164675. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCAGCCGACCGCTGCCACTATGGACTGCAACCAAGGAAAGTCTTCCCCAGTCCAAGGAGCAGTCGT
GTCTGACCTACATTGGGCTTTTCTCAGAACTTTGAACGATCCCATGCAAAGAATTCCCACCCTGAGGTG
GGTTACATACCTGCCAATGGCCAAAGGAACCGCCTTGTGAGACCAAATTGCTGACCTGGGTCAGTGCA
TGTGCTTTATGGTGTGGTGCATCTTTGGAGATCATCACCATATTTTACTTTTGGAGTCTTTAAAGAGG
AAGGGGAGTGGAGGGAACCCTGAGCTAGGCTTCAGGAGGCCCGCTCCTACGCAGGCTCTGCCACAGGG
GTTAGACCCAGGTCCGACGCTTACCTTCTGAGCTTCAAGTGCCTCCCCTATCAAATGAAGGGATG
GACAGCATGACCTCTGGGTGTCTCTCCAACCTACCAGTTCTAAAAAGGTATCAGATTCTATTGTGACT
TCATAGTGAGAATTTATTATAGATTATTTTTAGCTATTTTTCCATGTGTGAACCTTGAGTGATACTA
ATCATGTAAGTAAGAGTTCTCTTATGTATTATTTTCGGAAGAGGGGTGTGGTGACTCCTTTATATTCG
TACTGCACTTTGTTTTCCAAGGAAATCAGTGTCTTTACGTTGTTATGATGAATCCCACATGGGGCCG
GTGATGGTATGCTGAAGTTCAGCCGTTGAACACATAGGAATGTCTGTGGGGTGACTCTACTGTGCTTTA
TCTTTTAAACATTAAGTGCTTTGGTTTCAGAGGGGCAAGTATAAGCTCTGTTTCCCCTCTCCCCAAGC
CTTCAGCGAACGTGAAATGTGCGCTAAACGGGAAACCTGTTAATTCTAGATATAGGGAAAAAGGAAC
GAGGACCTTGAATGAGCTATATTCAGGGTATCCGGTATTTTGAATAGGGAAATAGGAAACCTTGTGGC
TGTGGAATATCCGATGCTTTGAATCATGCACTGTGTTGAATAACGTATCTGCTAAATCAGG
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001164675.2
Summary:	This gene encodes an enzyme that catalyzes the hydrolysis of sulfate esters by oxidizing a cysteine residue in the substrate sulfatase to an active site 3-oxoalanine residue, which is also known as C-alpha-formylglycine. Mutations in this gene cause multiple sulfatase deficiency, a lysosomal storage disorder. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2009]
Locus ID:	285362
MW:	37.8