

Product datasheet for SC211860

SC65 (P3H4) (NM_006455) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SC65 (P3H4) (NM_006455) Human 3' UTR Clone
Symbol:	SC65
Synonyms:	LEPREL4; NO55; NOL55; SC65
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006455
Insert Size:	1034 bp
Insert Sequence:	>SC211860 3'UTR clone of NM_006455 The sequence shown below is from the reference sequence of NM_006455. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GAGGCTGAGCCAGAGCCTGAACCTCGCATGAAGAGGGACACCCACACCGCTCAAGCTTGGGAAGCCTG
GTGCCGATGGCCCCACCCTCACCAGCCTGGGCAGCAGCAAGAATTTATTAATAAACTTAAGATGGGC
CAGGTGCGGTGGCTCACACCTGTAATCCCAGCATTTTGGGAGGCCAAGGTGGGTGGATCACTTGAGGCC
AGGAGTTCAAGACCAGCCTGGCCAACATGATGAGACCTCCGTCTCTACTAAAATACATAAATTAGCCGG
GTGTGGTGGCAGGCGCCTGAAATCCCAGCTACTCAAGAGGCTGAGGCAGGAGAATCGCTTGAACCTGGG
AGGCAAAGGTTGCAGTGAACCTGAGATTGCGCCACCGCACTCCAGCCTGGGCGACAGAGCGAGACTCCAT
CTTTAAAAAAAACAAGACGGGCGCGCACGGTGGCTCACGCCTGTAATCCCAGCACTGAGAGGCCGATC
ACTTGAGGTGAGGAGTTCAAGACCAGCCTGGCCAACATGGTGAAACCCCATCTCTACTAAAAATACAA
AAATTAGCCAGGCATGGTGGCACACCTGTAATCGTAGCTGAGGCAGGAGAATCGCCTGAACCCAGGA
GGCGGAGCTTGCAGTGAGCCGAGATCGTGCCACTGCACTCCAGCCTGGGCGACAGAGTGAGACTCCATC
TAAAAAAAAAAAAAAAAAACTTAAGATGGACACAGCTGACTGGACCCCATCCTGCCTCACCCATGGGT
GCTGCACCCAGACCCATCCTGCCACTTCTATGTCTCTGGACCACAGGATGGTGGTGGCATTGCAGGTT
GGCAAGTGGGCTGATGGGGTCCGCCCTCCTCACTGCTGAGCTCCTCACCTGGACAGTCTCCTGGACAAG
GAGTTTCCAGCTGCTGGGTGGAGTCTCAGGCCAAATTGCAGAGGGTCTCCAGGGTCTGAAGAGCACT
GGACTAAGAGTCTAGTGGTTCCAGGGCCCTGACCAGTAGGTGCTCAATAAATGTTTGTGTTGAATGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_006455.3</u>
Summary:	This nucleolar protein was first characterized because it was an autoantigen in cases on interstitial cystitis. The protein, with a predicted molecular weight of 50 kDa, appears to be localized in the particulate compartment of the interphase nucleolus, with a distribution distinct from that of nucleolar protein B23. During mitosis it is associated with chromosomes. [provided by RefSeq, Jul 2008]
Locus ID:	10609
MW:	37.9