

Product datasheet for **SC211873**

SLAIN1 (NM_001040153) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLAIN1 (NM_001040153) Human 3' UTR Clone
Symbol:	SLAIN1
Synonyms:	C13orf32
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001040153
Insert Size:	1039 bp
Insert Sequence:	>SC211873 3'UTR clone of NM_001040153 The sequence shown below is from the reference sequence of NM_001040153. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGC
GATGGAATTGGAGAGATGGTTGCTACTAATGCAGTTTTATGTACCCTTGAAAAATGGGAAAGAAGTAA
AAATGAGGGTTGTGTACCTAGCTGGCTGGGTAGCAGTGGATGTTGGGATATTCTTCCCTTTGTGTT
TTAATATATTTACTGCATTGTTTCTCAATGGACCAGTCACCAGAGACTAATTATGCACTTAAATATTT
GCCTGAGATACTGCAACATTCTCAAACCATGGTTGCAGTATTGTGACACTTAGATCTAGGAAGTTTTT
GTAGAAGTGTCTGTACCTGAATACTTTTTGAGAGAAATTAAGATGTATCAATAATGCTTTGCCATATGA
GTTTTTAAAGTAACTTGTTCATTTACTCACGTGTTCTAAACATCTTTCCATTACATGTTCTGTATTT
TAATACATTGCATATTGACAACCTAGTTCTATAATGTATGCTTTGAAATTTACTTTTTATAGTTTACA
GGAATTTTATTTTTGTGCCTATTTCTTTTACACCTATGTGAACCACTATGGAACAACCTAAATTTTG
TGCCATAAAAAATTTTTGTGGTAAGGTACTATTTTTTAGCTCTAGGGATATATCAGCAAAAACACAT
CATGCAATTTGAGACACATAATTTGTGTTGAATGAGCACAACATAATTTGAAGCATTGCAAGGAGATA
ACCAGACAGCAGAATTAATGGTCCTGTCTTTTCATTTTAAATTTATTTGTCATACATGGGTTTCATAT
TTATAACGGCATCATGAGCTCATTGCACTTAATACCTGCAATGTTGCTACTGTACCACAATTGATTTT
CAATACTTTATTACGAAGGATGAACTGTAATGTTTTTAAACAATGCTTCTGGAAATGAATGCATTTT
AAAGCAAATAAATCTTTTGATAGACCTTTTACAAAATCCATTTGCACTAATGAATGCTTTCTTATGGC
ATATACTTAATATTTGTACTGTGTACACTGCTGTTTTGGAATGTTGAGAAATAAGACTCTATTTCA
GCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001040153.4
Summary:	Microtubule plus-end tracking protein that might be involved in the regulation of cytoplasmic microtubule dynamics, microtubule organization and microtubule elongation. [UniProtKB/Swiss-Prot Function]
Locus ID:	122060
MW:	40.9